


# Fc Capture A/G Kit

with red dye **Ra**

Dynamic Biosensors GmbH & Inc.  
HK-FC-3 v4.1




## Key Features

- This kit is designed for capture of **IgG molecules** binding to **Protein A and G**.
- Compatible with **helix<sup>®</sup> Adapter Chip**.
- Includes **Adapter strands** and **Ligand strand** modified with Protein A/G for **20 regenerations**.
- For functionalization of **Spot 1** and **Spot 2**.
- **Adapter strands 1** and **2** carry a moderately hydrophilic red dye (**Ra**) with a single positive net charge.

## helix<sup>®</sup> Adapter Chip Overview

2 spots with 2 different anchor sequences for DNA-encoded addressing. Spot 1 is functionalized with the capture molecule while Spot 2 is used as real-time reference.



## Product Description

Order Number: **HK-FC-3**

Table 1. Contents and Storage Information


Material	Cap	Concentration	Amount	Buffer	Storage
<b>Protein A/G - Ligand strand</b>	Red	500 nM	2 x 100 µL	TE40 <sup>[1]</sup>	2-8°C
<b>Adapter strand 1 - Ra</b>	Black	400 nM	2 x 100 µL	TE40 <sup>[1]</sup>	-20°C
<b>Adapter strand 2 - Ra - lfs</b>	White	200/250 nM	2 x 200 µL	TE40 <sup>[1]</sup>	-20°C

For research use only.

This product has a limited shelf life, please see expiry date on label.  
After preparation of ready to use solution the expiry date is **6 months**.

Protein A/G binds to all IgG subclasses from various mammalian species, including all IgGs that bind to both Protein A and Protein G.

## Workflow of a heliX<sup>®</sup> FC - capture assay



1. The anchor strand (ssDNA) immobilized to the surface of the **heliX<sup>®</sup> Adapter Chip** is hybridized with complementary DNA strands modified with **Protein A/G**.
2. The **IgG** of interest is captured on the surface during the measurement run.
3. Measurement of the analyte binding kinetics.
4. Surface regeneration by injection of a high pH solution. Chip surface goes back to the original state. This step can be followed by a new hybridization of fresh ligand with Protein A/G.

## Preparation

1. Mix 100  $\mu$ L **Protein A/G - Ligand strand** with 100  $\mu$ L **Adapter strand 1 - Ra**.
2. Incubate the solution of step 1 at **RT** at **600 rpm** for **30 min** to ensure complete hybridization.
3. Mix 200  $\mu$ L **Adapter strand 2 - Ra - lfs** to the sample after step 2.

The solution (400  $\mu$ L in total) is ready to use for a biochip functionalization.

Please aliquot and store the ready to use solution at 2-8°C. **Use up within 6 months.**

The kit contains material for the preparation of two separate ready to use solutions with 400  $\mu$ L each.

## Assay Setup in heliOS

Go to **heliOS** > create a **New Assay Workflow** > add **Custom Assay** > load **Capture with Kinetics** > modify the parameters based on your needs and run the assay.

Suggested assay parameters (e.g., flow rate, time, LED power, etc.) are within the **heliOS** assay.

### TIP

*If the fluorescent change signal upon analyte binding is very small, consider using the conjugation approach, in order to move the binding site closer to the dye.*

For further questions, please contact the support team at [support@dynamic-biosensors.com](mailto:support@dynamic-biosensors.com)

## Useful Order Numbers

Table 2. Order Numbers

Product Name	Comment	Order No
<b>heliX<sup>®</sup> Adapter Chip</b>	Chip with 2 detection spots	ADP-48-2-0
<b>10x Passivation solution</b>	For passivation of chip surface	SOL-PAS-1-5
<b>Regeneration solution</b>	For regeneration of chip surface	SOL-REG-1-5

## Contact

**Dynamic Biosensors GmbH**

Perchtinger Str. 8/10  
81379 Munich  
Germany

**Dynamic Biosensors, Inc.**

300 Trade Center, Suite 1400  
Woburn, MA 01801  
USA

**Order Information** [order@dynamic-biosensors.com](mailto:order@dynamic-biosensors.com)

**Technical Support** [support@dynamic-biosensors.com](mailto:support@dynamic-biosensors.com)

[www.dynamic-biosensors.com](http://www.dynamic-biosensors.com)

Instruments and chips are engineered and manufactured in Germany.

©2024 Dynamic Biosensors GmbH | Dynamic Biosensors, Inc. All rights reserved.

[1] TE40: 10 mM Tris, 40 mM NaCl, 0.05 % Tween20, 50 µM EDTA, 50 µM EGTA