

## Normalization solution (green dye)

for RT-IC measurements in the green channel

Dynamic Biosensors GmbH & Inc.

NOR-Ga v1.0



## Key Features

- For the normalization of the fluorescent signals on Spot 1 and Spot 2 of a **heliX<sup>cyto</sup>** chip
- Enables correct real-time referencing of the **green fluorescent signals** during **RT-IC** measurements
- Compatible with all **heliX<sup>cyto</sup>** chips
- The **Normalization solution (green dye)** contains a hydrophilic green dye with a single negative net charge

## Product Description

Order Number: **NOR-Ga**

*Table 1. Contents and Storage Information*

Material	Cap	Concentration	Amount	Storage
<b>Normalization solution-Ga</b>	Green	10 $\mu$ M	6x 100 $\mu$ L	-20 $^{\circ}$ C

For research use only.

This product has a limited shelf life, please see expiry date on label.

## Preparation

Use this green dye normalization solution for RT-IC measurements in the green channel (analyte label-dependent).

**Dilute** the 10  $\mu$ M normalization stock solution to a working concentration with **running buffer**.

The concentration of the normalization solution should approximately correspond to the fluorophore concentration in the highest analyte concentration to be measured. This can be calculated using following equation:

$$c_n[M] = c_f[M] = c_a \cdot DOL$$

$c_n$ : Concentration of the normalization solution in the desired color

$c_f$ : Concentration of dye in the labeled analyte solution

$c_a$ : Highest concentration of analyte that should be measured

$DOL$ : Degree of labeling (ratio of dye to analyte)

Diluted solutions can be stored at 2-8°C for up to 7 days.

## Application Note

In the RT-IC measurement, the fluorescent signal of the normalization solution should be in a **similar range** as the highest signal coming from bound analyte (raw data).

The absolute fluorescent signal is dependent on normalization solution concentration and the excitation power applied in the measurement. The excitation power has to be selected based on the following parameter:

a. Fluorophore concentration in analyte solution:

The fluorophore concentration depends on the analyte concentration used in the measurement as well as the degree of labeling of the analyte. For high DOL and high analyte concentrations, lowering the excitation power might be required.

b. Expected binding signal:

Highly expressed targets on a cell can bind more molecules of labeled analyte. In case of highly overexpressed targets a strong binding signal can be expected. To avoid the shutter closing, lowering the excitation power might be considered.

c. Chip type:

Different chip types have varying fluorescent background. The bigger the traps and the more traps on the chip, the higher the background signal. Therefore, L5 chips might require lower excitation power than applied to M5 chips.

For a **starting point** of excitation power and norm. solution concentration to be used in an RT-IC experiment, please refer to Table 2.

Table 2. Relation of fluorophore concentration, normalization solution concentration, and excitation power suitable for a heliX<sup>cyto</sup> M5 chip

Analyte dye conc. = analyte conc x DOL	Excitation power	Concentration Normalization solution	Dilution Normalization solution
25 nM	0.5	25 nM	1:400
50 nM	0.3	50 nM	1:200
100 nM	0.2	100 nM	1:100
300 nM	0.1	300 nM	1:33
500 nM	0.08	500 nM	1:20
1 µM	0.05	1 µM	1:10
2.5 µM	0.02	2.5 µM	1:4

Note: This table is for your guidance. However, the final signal recorded in the heliX<sup>cyto</sup> depends on many factors. Thus, some optimization will be required for each system.

## Contact

**Dynamic Biosensors GmbH**

Perchtinger Str. 8/10  
81379 Munich  
Germany

**Dynamic Biosensors, Inc.**

300 Trade Center, Suite 1400  
Woburn, MA 01801  
USA

**Order Information** [order@dynamic-biosensors.com](mailto:order@dynamic-biosensors.com)

**Technical Support** [support@dynamic-biosensors.com](mailto:support@dynamic-biosensors.com)

[www.dynamic-biosensors.com](http://www.dynamic-biosensors.com)

Instruments and chips are engineered and manufactured in Germany.

©2024 Dynamic Biosensors GmbH | Dynamic Biosensors, Inc. All rights reserved.