

## Real-time measurement of Ni<sup>2+</sup>-ion chelation kinetics

Key words: NTA, His-tagged proteins, Chelation, Ions


### Background

Nitrilotriacetic acid (NTA) is a common chelating agent, which forms coordination compounds with metal ions (chelates) such as Ni<sup>2+</sup>, Cu<sup>2+</sup> and Co<sup>2+</sup>. It is primarily used for the purification and immobilization of His-tagged proteins, which are captured by NTA-Ni<sup>2+</sup> chelates bound to solid supports. The loading of NTA with Ni<sup>2+</sup>-ions is an essential step for the activity of the chelate. Thus, we determined the association rate of Ni<sup>2+</sup> to NTA and investigated the stability of the chelate against a competitive chelator, EDTA. Such information is required for research applications involving His-tag capture and industrial protein purification methodologies.

### Methods and Results

Using a standard switchSENSE chip with 48 bp DNA, the binding of Ni<sup>2+</sup> to NTA was measured by Fluorescence Proximity Sensing. Complementary trisNTA-modified DNA from the Dynamic Biosensors NTA kit was used to functionalize the surface with NTA groups.


When chelated by the NTA groups, the bound Ni<sup>2+</sup>-ions quench the dye fluorescence by approximately -15% due to the Fluorescence Proximity effect. Binding kinetics follow a single-



**Figure 1** | Real-time measurement of binding of Ni<sup>2+</sup>-ions to NTA measured by Fluorescence Proximity Sensing. Red lines are single exponential fits with an association rate of  $k_{ON} = 1.2 \pm 0.1 \times 10^4 \text{ M}^{-1}\text{s}^{-1}$ .

exponential behavior (Fig. 1) with an association rate of  $k_{ON} = 1.2 \pm 0.1 \times 10^4 \text{ M}^{-1}\text{s}^{-1}$ .

Ni<sup>2+</sup>-ions are stably bound to NTA for many hours. To remove Ni<sup>2+</sup> from NTA, high concentrations of a competitive chelator like EDTA are necessary. Only when using EDTA concentrations higher than 1 mM the dissociation of Ni<sup>2+</sup> could be observed (Fig. 2) and even at an EDTA concentration of 10 mM the dissociation rate was relatively low:  $k_{OFF} = 6.2 \pm 0.4 \times 10^{-3} \text{ s}^{-1}$ .



**Figure 2 |** Real-time measurement of the dissociation of Ni<sup>2+</sup>-ions from NTA via competition with EDTA. The red line is a single exponential fit. The dissociation rate at 10 mM EDTA is  $k_{OFF} = 6.2 \pm 0.4 \times 10^{-3} \text{ s}^{-1}$ .

## Conclusions

The binding kinetics of Ni<sup>2+</sup>-ions to NTA could be measured using the Fluorescence Proximity Sensing measurement mode. For complete saturation of an NTA-surface within 1 minute, a concentration of at least 10  $\mu\text{M}$  Ni<sup>2+</sup> is required. A low background concentration of e.g. 100  $\mu\text{M}$  EDTA does not affect the chelated ions. To wash Ni<sup>2+</sup> from NTA, EDTA concentrations of at least 10 mM must be used.

Author | Andreas Langer

Contact | [info@dynamic-biosensors.com](mailto:info@dynamic-biosensors.com)

Dynamic Biosensors GmbH | January 2016