


Biophysical Characterization of Liposomes and Antibody-Receptor Interactions on VLPs with the switchSENSE® Biosensor


Hanna Mueller-Landau, Ulrich Rant and Wolfgang Kaiser; Dynamic Biosensors GmbH, Martinsried, Germany

Virus-Like-Particles (VLPs) are membrane bound vesicles that are an attractive vaccine platform authentically resembling their host virus in structure and antigen expression. We report a novel automated methodology for characterizing VLPs with electro-switchable biosurfaces (switchSENSE®). Stable surface functionalization with VLPs (200 nm in diameter) using antibody modified DNA monolayers enables characterization of molecular interactions towards membrane proteins in physiological conditions. Additionally, a strategy to immobilize liposomes on the sensor surface using cholesterol modification is presented. Temperature dependent phase transitions of immobilized liposomes was monitored in real-time. Information about specific lipid melting temperatures is crucial for the application of liposomes as thermosensitive drug-delivery systems.

Strategies for Surface Functionalization with Virus-Like Particles




Real-Time Binding Analysis of Antibodies against Membrane Receptors



Detecting Antibody-Receptor Interactions with switchSENSE®

- ✓ Stable surface immobilization of VLPs via DNA encoded addressing
- ✓ Independent monitoring of multi-step assays on one detection spot using the DRX² dual color system
- ✓ Time-resolved observation of antibody binding towards membrane receptors

Characterization of Liposomes on Electro-Switchable Biosurfaces




Characterizing Liposomes with switchSENSE®

- ✓ Stable surface immobilization of liposomes via DNA encoded addressing
- ✓ Time-resolved analysis of liposome binding and lipid phase transitions
- ✓ Characterization of specific melting temperatures of lipid formations

Electro-Switchable DNA Nanolevers


switchSENSE® | Measurement Cycle

Regeneration of sensor surface by DNA denaturation and repeated hybridization



switchSENSE® | Measurement Modes


Two independent measurement modes for comprehensive signal acquisition



switchSENSE® | DRX² Instrument

DRX² instrument featuring two light sources & two photon counters optimized for red & green fluorophores for dual binding analysis.

Multi-electrode biochip with 4 separate flow channels each comprising 6 detection spots.



- High Sensitivity**
Measurement of analyte concentrations from fM to mM.
From ultra-fast to ultra-slow kinetics.
LOD = 10 fM
- Kinetics and Affinity**
Determination of binding rate constants k_{ON} , k_{OFF} and dissociation constants K_D in real-time.
- Size and Conformation**
Analysis of protein diameters on chip with 0.1 nm accuracy, and monitoring of conformational changes.
- Cooperativity and Avidity**
Identification of multiple binding sites in a single measurement using variable capture molecule densities.
- Thermodynamics**
Characterization of melting transitions, thermal stability or thermodynamic analyses. $8^\circ - 75^\circ\text{C}$