

10x Buffer PE140 pH 7.4

After dilution this buffer is intended for use as a running buffer for the **heliX**[®] and DRX series.

Product Description

Order Number **BU-PE-140-10**

Material	Amount	Storage	Comment
10x Buffer PE140 pH 7.4 (100 mM Na ₂ HPO ₄ /NaH ₂ PO ₄ , 1.4 M NaCl, 500 µM EDTA, 500 µM EGTA and 0.5 % Tween20; 0.2 µm sterile filtered)	50 mL	r.t.	

For *in vitro* use only.

This product has a limited stability, please see expiry date on label.

At temperatures lower than 18°C precipitation may occur.

Preparation

Dilute 10x Buffer PE140 pH 7.4 with deionized, 0.2 µm filtered water by mixing the complete solution (50 mL) with 450 mL ultrapure water. After dilution PE140 Buffer is ready for use (10 mM Na₂HPO₄/NaH₂PO₄ pH 7.4, 140 mM NaCl, 50 µM EDTA, 50 µM EGTA and 0.05 % Tween20).

The diluted buffer should be stored at 2-8°C.

Contact

Dynamic Biosensors GmbH

Perchtinger Str. 8/10
81379 Munich
Germany

Phone: +49 89 89 74 544 0

Dynamic Biosensors Inc.

300 Trade Center, Suite 1400
Woburn, MA 01801
USA

Phone: +1 781 404 6126

Order Information


Phone: +49 89 89 74 544 0

Email: order@dynamic-biosensors.com

Technical Support

Phone: +49 89 89 74 544 66

Email: support@dynamic-biosensors.com



switchSENSE® is a proprietary measurement technology by Dynamic Biosensors GmbH. Instruments and biochips are engineered and manufactured in Germany.

©2023 Dynamic Biosensors GmbH | Dynamic Biosensors Inc. All rights reserved.