


switchSENSE® Demo Kit

DK-Kinase: Protein – Small Molecule (k_{ON} | k_{OFF} | K_D)

Association and Dissociation of a Small Molecule

Aims

- Familiarization with the software tools for performing switchSENSE® experiments
- Determination of the association and dissociation rates of a small molecule (k_{on} , k_{off} , K_D)



Workflow

Design your experiment



switchBUILD

Step-by-step guidance through the workflow

Perform your experiment



switchCONTROL

Finetuning of the experiment (if needed)

Analyze your results



switchANALYSIS

Simple drag'n'drop format for efficient data interpretation

Product description

Product Code DK-KS-1-B96

Measurement Time ~2 h 04 min

In this experiment, the binding of the small molecule Staurosporine to the Kinase PKA is analyzed. We start with immobilizing the Y-structure on the biosensor surface. Next, the association and dissociation of three consecutive concentrations of Staurosporine to PKA is observed.

Before getting started, transfer one aliquot of each sample listed in Table 1 to the fridge to ensure slow thawing.

TABLE 1 | Contents and storage information



Material	Cap	Amount	Storage	Comments	
cNL-A96 / cNL-B96 (500 nM each)	blue	3 x 40 µL	-20°C	in PE40	
Y-structure-PKA (200 nM)	green	3 x 40 µL	-20°C	in TE40, 0.5 mM TCEP	
Staurosporine	red	5 nM	3 x 500 µL	-20°C	in Running Buffer
		10 nM			
		20 nM			


TABLE 2 | Additional material required

Material	Cap	Amount	Storage	Comments
MPC-96-2 biochip (red amplitude > 400 kcps)	-	1	2-8°C	use in DRX or DRX ²
Kinase Buffer (1x): 10 mM TRIS-HCl, pH 7.4, 40 mM NaCl, 2 mM MgCl ₂ , 1 % DMSO, 50 µM EDTA, 50 µM EGTA, 0.05 % Tween 20, sterile filtered	transparent	2 x 50 mL	2-8°C	Running buffer
TE40 buffer (1x): 10 mM TRIS-HCl, pH 7.4, 40 mM NaCl, 50 µM EDTA, 50 µM EGTA, 0.05 % Tween 20, sterile filtered	transparent	2 x 50 mL	2-8°C	Aux buffer
DI Water	transparent	1 x 12 mL	2-8°C	
Passivation solution (10x)	yellow	1 x 1 mL	2-8°C	dilute in 9 ml DI water
1x Regeneration solution	transparent	1 x 12 mL	2-8°C	


Large glass vials (10 mL), small glass vials (1.5 mL), caps required. For *in vitro* use only.

Planning the experiment

1. Start switchBUILD by clicking on the  icon on the desktop and click on **+ New** to create a new assay. The interface is arranged in two different columns: *Assay* and *Properties*. The *Assay* column shows all elements of your assay. By default, it includes a panel to choose the basic settings of your assay and a *Passivation* element, which should be run before the measurement to prevent unspecific binding to the biosensor surface. Next to *Assay* the estimated runtime of your assay is shown. In the *Properties* column, the individual measurement parameters for each assay element can be set. The info button  always provides further information about the current topic.
2. Now, let's start with defining the general *Properties* of your assay by selecting the first panel of your assay. Choose the **flow channel** you want to run the experiment in by clicking on the channel – the selected channel is highlighted in green. Furthermore, the running buffer *Buffer* and auxiliary buffer *Aux* can be specified. For this demo select “**Custom**” as *Buffer*. You can enter your buffer composition if wanted. To specify the type of Biochip used, select “**MPC-96-2**” from the *Select chip* list, which indicates that it is a Multi Purpose Chip (MPC) with a functionalized DNA-length of 96 base pairs and two different DNA sequences. The two sequences are spatially separated with sequence A on electrodes 1 and 2, and sequence B on electrodes 3-6. You can use within one measurement sequence A as control and sequence B for the specific sample of interest.



3. To add assay elements to the sequence, click  in the *Assay* column and choose the desired element. As the goal of this experiment is to determine binding kinetics, we need a **Kinetics** element.



4. In the next step, you should define the *Properties* of the Kinetics assay. In the upper section general information about your interaction experiment is requested. Select **Conjugate Hybridization** for *Immobilization*, which assumes that your ligand is conjugated to the complementary DNA. Change the *Measurement* type to **Static Mode**. In this mode, the DNA is kept in a defined upright position. Update the name of your immobilized *Ligand* to **Y-Structure-PKA**, since we are hybridizing the Protein Kinase A (PKA) conjugated to DNA oligos forming a Y-Structure, and select **200 nM** as concentration. The molecular weight of the kinase is around **40 kDa**. The name of your *Analyte* is **Staurosporine**; the small molecule has a molecular weight of **0.5 kDa**.

5. In the *Interaction* tab you are asked to give an estimation of the expected interaction parameters (k_{on} , k_{off} , K_D). You can either choose saved settings via **PRESETS** or enter custom values. Please enter the following parameters, either by moving the sliders or by typing the values in the corresponding fields: $K_D = 2.0E-9$ M, $k_{on} = 7.0E+06$ M⁻¹s⁻¹, $k_{off} = 1.4E-02$ s⁻¹. When using the sliders, one of the three values will always be fixed, as each value is dependent on the other two: $K_D = \frac{k_{off}}{k_{on}}$. The values entered are used to predict the Fraction Bound [%] shown in the graph on the right, which is updated on the fly when you make changes in the next step:
6. As *Experimental Parameters* we keep **3 concentrations** of Staurosporine, starting from the stock concentration of **2E-8 M**, which will be diluted by a *factor* of **2**. The resulting concentrations are displayed below the graph on the right. Additionally, add buffer control runs by activating the checkboxes next to **with blank run** for both Association and Dissociation. Next enter the *Association volume* of **460 μL**, with a pump rate of **500 μL/min**. The time is the result of volume divided by pump rate and can be used for setting the volume as well. The *Dissociation* should be run with **12000 μl** at a pump rate of **1000 μl/min** with the peristaltic pump (PP). If there is only a syringe pump (SP) available, you can use 2400 μL with a flowrate of 200 μL/min for dissociation.

7. In the third section, you are asked to select the measurement spots that will be used in your experiment. The selected spots are colored in green. For this demo, choose **one measurement spot** out of spots 3-6 and select it for both *Association* and *Dissociation*. **Deselect** all electrodes for the *Stopped Flow*.

In addition, you can decide if the flow channel should be functionalized with fresh ligand after each concentration of analyte. This includes a regeneration, during which the biochip is washed with a proprietary regeneration solution (SOL-REG-12-1) to ensure a fully single-stranded state of all covalently immobilized DNAs. Subsequently, fresh conjugate DNA is re-hybridized. It is advisable to perform regenerations between concentrations if no complete dissociation of the analyte has taken place beforehand. Since Staurosporine is dissociating within 10 - 15 min, instead of regenerating we will fully dissociate each analyte concentration. To do so, please **uncheck** the two checkboxes *with Regenerations* and *only last concentration*. To improve the signal-to-noise ratio, you can increase the excitation power by clicking on **ADVANCED** and selecting the excitation power **2 LED** for *Association* and *Dissociation*.

Measurement Spots

Association	1	2	3	4	5	6
Stopped Flow	1	2	3	4	5	6
Dissociation	1	2	3	4	5	6

with Regenerations

use Dynamic Mode

only last concentration

ADVANCED

- At last, add a **Standby** assay element to your *Assay* sequence and choose **Reset surface** only. This step removes the ligand, functionalizes the chip with bare DNA, fills the measurement channel with air and thereby prepares the biochip for storage after your measurement.




- Save your newly designed assay as a DBB file by clicking on and name the file "DK-Kinase_1".


- You will be taken to the *Autosampler* tab, where you can see how the autosampler of the DRX should be loaded for the experiment. In the middle column you can find a list of all solutions and their concentrations needed for your measurements. When you click on an item in the solution list, its respective position in the autosampler display (right column) is indicated with a green circle.

Sample	Volume (µl)	Conc (M)	Type	Vial
A1 Waste	-440		Waste	Large
B1 Water	filled		Water	Large
C1 Passivation solution	1610		Solution	Large
A2 X40	1200		Solution	Large
A5 Regeneration solution	65		Solution	Small
B5 cNL-A06 and B96	40	5.00E-07	Ligand	Small
C5 Y-structure-PKA	40	2.00E-07	Ligand	Small
A10 Staurosporine	500	0	Blank	Small
B10 Staurosporine	500	5.00E-09	Analyte	Small
C10 Staurosporine	500	1.00E-08	Analyte	Small
D10 Staurosporine	500	2.00E-08	Analyte	Small


Congratulations, you just designed your own **switchSENSE®** experiment!

Conducting the experiment

1. Start **switchCONTROL** by clicking on the  icon on the desktop and enter **your name** when prompted (the measured data will be automatically saved in a folder with your name).
2. To insert a new chip, click on **Chip** → **Eject Chip Holder** in the upper left corner of the software window to eject the chip holder from the DRX.



3. The following message will pop up. Indicate which channels have been used on the inserted biochip by checking their respective checkboxes, and then click **Ok**.






4. After venting, the chip holder will be ejected. Open the chip holder by pushing both white sliders on the right to the back. Use the vacuum pen to remove the chip currently in place and insert your **MPC-96-2** chip. Make sure to place it with the label facing upwards. Once the chip is in place, close the chip holder while pushing the white slider back.



5. Click on **Chip** → **Retract Chip Holder** to retract the new biochip into the DRX and follow the instructions of the **Retract assistant**. Enter the Chip-ID as shown on the packaging of your biochip. A message will pop up asking if the chip should be **aligned** – the alignment is necessary to precisely position the electrodes for measurement. An alignment should be performed every time a chip is placed into the DRX (except for a Cleaning chip which is not used for measurements). The next message will ask you if you wish to **prime** the microfluidics system with the running buffer – the priming should be done if this is your first measurement of the day or if your running buffer is different from that of the prior measurement. Make sure to connect the *Running Buffer* (Kinase Buffer) to the **Buffer** and the *Aux Buffer* (TE40 buffer) to the **AUX** reservoirs on the top left (exterior) of the instrument. Finally, you are asked if you want to check the status of the electrodes and regenerate the channel before starting the experiment. These steps are not necessary for the purposes of this demo and the default settings can be maintained. The alignment and priming will take a few minutes. You can follow the alignment by watching the small camera window on the top left of the software window, where you will see the focusing on the four alignment crosses.



- While the alignment and priming are taking place, please prepare the samples as indicated in switchBUILD.
- Load your assay *DK-Kinase_1*** into switchCONTROL by clicking on  at the left of the software window {1}. Next, click on the **Autosampler** tab {2} to view the sample positions listed within the script, and load the sample solutions into the autosampler accordingly. Close the autosampler door and start the measurement by clicking on  {3}.




The screenshot displays the software interface with three main sections:



- Taskflow:** A list of tasks including C1, Passivation, Kinetics, Standby, and Deinitialize. A folder icon is highlighted with a blue box and the number 1.
- Autosampler:** A grid representing the sample plate layout. The grid is labeled with columns A-F and rows 1-10. A blue box with the number 2 highlights the Autosampler tab.
- Autosampler Sample Preparation:** A table with columns: No., Sample, Volume (µl), Conc (M), Vial, Inset, and Color. The table contains data for rows A10, B10, C10, and D10, all with 'Staurosporine' as the sample and a volume of 500 µl. A 'Start/Stop' button is highlighted with a blue box and the number 3.

No.	Sample	Volume (µl)	Conc (M)	Vial	Inset	Color
A 10	<input checked="" type="checkbox"/> Staurosporine	500	0.00E+0	small		Orange
B 10	<input checked="" type="checkbox"/> Staurosporine	500	5.00E-9	small		Orange
C 10	<input checked="" type="checkbox"/> Staurosporine	500	1.00E-8	small		Orange
D 10	<input checked="" type="checkbox"/> Staurosporine	500	2.00E-8	small		Orange
E 10	<input type="checkbox"/>	0				Black
F 10	<input type="checkbox"/>	0				Black
A 9	<input type="checkbox"/>	0				Black
B 9	<input type="checkbox"/>	0				Black
C 9	<input type="checkbox"/>	0				Black
D 9	<input type="checkbox"/>	0				Black
E 9	<input type="checkbox"/>	0				Black
F 9	<input type="checkbox"/>	0				Black

Measurement Steps

If working with a DRX² DualColor instrument, choose to display the **red dye only** by clicking on the color block next to *Display* above the measurement plots .

- Passivation and Status


Before starting an experiment, the biosensor surface needs to be passivated to prevent unspecific surface binding, and thus enable maximal DNA-switchability. During the passivation step, the chip is incubated in passivation solution and the DNA levers are switched with a frequency of 0.2 kHz. Positive potential leads to a low fluorescence signal (as the fluorophore lies close to the surface), while negative potential leads to a high fluorescence signal (as the fluorophore is positioned away from the surface). The *passivation* step takes 10 minutes and should be run completely if the channel had not been used on the same day. If the channel has already been used, you can skip this passivation step after a few minutes by clicking on . To do so, the taskflow needs to be unlocked first by clicking on .

After the passivation the chip will be subjected to a *status* diagnosis: Each electrode is calibrated (The change of fluorescence intensity in response to the applied voltage is recorded and the inflection point determined. More positive voltages attract the negatively charged DNA and orient the DNA nanolevers closer to the quenching gold surface.) and the DR and Fluorescence values are recorded. Use this information in case electrodes need to be changed.


Please make sure to use an electrode, for which the fluorescence amplitude is above 400 kcps.

- Functionalization

You will be able to follow the functionalization of the chip with the PKA-conjugated Y-structure in real-time in the **red channel**. Before the *functionalization*, the dsDNA is denatured by the injection of a *regeneration* solution (high pH), therefore the measurement starts with the switching of ssDNA. After a short while, a mix of the fresh oligos is injected into the flow channel and the Y-structure is formed. As dsDNA




is more rigid than ssDNA, the fluorophore is pushed away further from the surface during switching, which normally leads to a fluorescence increase. In the case of the Y-structure, a competitor compound is hybridized next to the dye (blue), which efficiently quenches the fluorophore (fluorescence decreases). Since this effect is more pronounced than the fluorescence increase caused by elongation of dsDNA, the overall fluorescence decreases during hybridization. The functionalization takes about 15 minutes for completion.




- Association and Dissociation of Staurosporine

The binding of Staurosporine to PKA is measured in Static Mode, which implies that the DNA nanolevers are not actuated but kept in an upright position. Since without switching of the DNA nanolevers no Dynamic Response can be determined, the plot for Dynamic Response is only a flat line. For the same reason, only the Fluorescence in the top position is plotted below.


For 30 secs, the fluorescence intensity of the nanolevers functionalized with PKA will be measured. After that, the Staurosporine solution is injected. Now you can follow the binding of the small molecule to PKA in real-time: The quenching effect of the immobilized compound (blue) is strongest, if it is not bound to the kinase (orange). That state is induced by the analyte, here Staurosporine (grey), when it associates to the kinase. Thus, a decrease in fluorescence can be observed until equilibrium is reached. After each concentration, running buffer is injected into the microfluidic channel. You can follow the dissociation of the small molecule in real-time by an increase in fluorescence. Since Staurosporine dissociates almost completely within a few minutes, no fresh conjugate is immobilized between the different concentrations of Staurosporine.






Association




Dissociation



Data analysis with the switchANALYSIS software

1. Start switchANALYSIS from the  icon on the desktop and click on  New to create a new assay. Now you can load your data by clicking on **IMPORT EXPERIMENT FILE** in the lower left corner and selecting your zipped measurement data, which is saved in > folder with your username > subfolder with the date.
2. You will find the assay structure of your switchBUILD script (Passivation – Kinetics – Standby) when expanding the measurement file (indicated by the start time of your measurement) by **clicking on the triangles**. Further expansion of the *Kinetics* dataset {1} reveals data for the Functionalization and the Association and Dissociation of Staurosporine. Drag-and-drop the folders **Association** and **Dissociation** into the right column to select the data {2}. Both the red and green channel recorded data during the measurements, however, we did not use any green-labeled DNA, therefore **delete all green** marked data with  {3}. You can also mark all green measurements in the right column by pressing *Ctrl* on the keyboard while selecting and delete them simultaneously with *Del*.







The screenshot shows the software interface with a table of experimental data and a list of measurements on the right. The table has columns for Name, Time, Ligand, Sample, Concentration, Electrode, Color, and Type. The 'Kinetics' folder is expanded, showing sub-folders for Functionalization, Association, and Dissociation. An orange box labeled '1' highlights the 'Kinetics' folder. An orange arrow labeled '2' points to the 'Association' and 'Dissociation' folders being dragged into the right column. The right column shows a list of measurements with checkboxes and color indicators (red and green). An orange box labeled '3' highlights a green measurement, indicating it should be deleted.


Name	Time	Ligand	Sample	Concentration	Electrode	Color	Type
2018-07-25 09:01-DK-Kinase_1	2018-07-25 09:04						
Passivation	2018-07-25 09:04	Air		0 M	●●●●●●●●	●●	
Kinetics	2018-07-25 09:31				○●○●○●○●	●●	
Functionalization	2018-07-25 09:31	NL-896	Y-structure	200 nM	○●○●○●○●	●●	
Association	2018-07-25 09:40	Y-structure	Staurosporine		○●○●○●○●	●●	
Dissociation	2018-07-25 09:43	Y-structure	Staurosporine		○●○●○●○●	●●	
Standby	2018-07-25 10:55	NL-896	cNL-A96 ar	500 nM	○●○●○●○●	●●	


3. Click on **CREATE NEW ANALYSIS** in the lower right corner and select the **Kinetics** analysis type to analyze the kinetic rates.

Please select an analysis type



- On the **SETTINGS** page you will see the data you just imported. Association and dissociation should be paired automatically in the main column, or you can drag-and-drop them from the left column {1} into the respective box and add them with the middle button Υ {2}. The fitting range is preselected automatically as well and indicated by the sliders. Select each measurement one by one in order to verify that the fitting borders are set accurately. They can be adjusted by entering a number after double-clicking on the value or by moving the sliders {3}. Proceed by clicking on **FIT ANALYSIS** in the lower right corner.



- You will see a mono-exponential fitting curve for the association and dissociation data. To incorporate all data except the buffer control, the global mono-exponential fitting model is selected in the left column. One value for k_{on} , k_{off} and K_D is calculated from all curves.
- Finally, you can save your analysis by clicking on .

Congratulations, you are a **switchSENSE[®]** expert user now!


Contact

Dynamic Biosensors GmbH
Lochhamerstr. 15
82152 Martinsried/Planegg
Germany

Dynamic Biosensors Inc.
9705 Carroll Centre Road, Suite 100
San Diego, CA 92126
USA

Order Information Phone: +49 89 89 74 544 0
Email: order@dynamic-biosensors.com

Technical Support Phone: +49 89 89 74 544 66
Email: support@dynamic-biosensors.com



switchSENSE® is a proprietary measurement technology by Dynamic Biosensors GmbH.

Instruments and biochips are engineered and manufactured in Germany.

© 2018 Dynamic Biosensors GmbH | Dynamic Biosensors Inc. All rights reserved.

Munich DE | San Diego, CA USA | Betchton UK | Paris FR | Tokyo JP | Singapore SG