


1x Buffer C pH 8.0

coupling buffer for nanolever conjugation

Dynamic Biosensors GmbH & Inc.

BU-C-150-1 v2.1



Product Description

Order Number: **BU-C-150-1**

Table 1. Contents and Storage Information

Material	Composition	Amount	Storage
1x Buffer C pH 8.0	50 mM Na ₂ HPO ₄ /NaH ₂ PO ₄ , 150 mM NaCl; 0.2 µm sterile filtered	50 mL	2-8°C

For research use only.

This product has a limited shelf life, please see expiry date on label.

Contact

Dynamic Biosensors GmbH

Perchtinger Str. 8/10
81379 Munich
Germany

Dynamic Biosensors, Inc.

300 Trade Center, Suite 1400
Woburn, MA 01801
USA

Order Information order@dynamic-biosensors.com

Technical Support support@dynamic-biosensors.com

www.dynamic-biosensors.com

Instruments and chips are engineered and manufactured in Germany.

©2024 Dynamic Biosensors GmbH | Dynamic Biosensors, Inc. All rights reserved.