


1X BUFFER C pH 8.0

coupling buffer for nanolever conjugation

Dynamic Biosensors GmbH
BU-C-150-1 v2.1



Product Description

Order Number: **BU-C-150-1**

Table 1. Contents and Storage Information

Material	Composition	Amount	Storage
1x Buffer C pH 8.0	50 mM Na ₂ HPO ₄ /NaH ₂ PO ₄ , 150 mM NaCl; 0.2 µm sterile filtered	50 mL	2-8°C

For research use only.

This product has a limited shelf life, please see expiry date on label.

Contact

Dynamic Biosensors GmbH

Perchtinger Str. 8/10
81379 Munich
Germany

Bruker Scientific LLC

40 Manning Road, Manning Park
Billerica, MA 01821
USA

Order Information order.biosensors@bruker.com

Technical Support support.dbs@bruker.com

www.dynamic-biosensors.com

Instruments and chips are engineered and manufactured in Germany.

©2025 Dynamic Biosensors GmbH

For Research Use Only. Not for use in clinical diagnostic procedures.