

## Cleaning Solution 1

Detergent-based solution for everyday cleaning procedures of heliXcyto and chips.

Dynamic Biosensors GmbH & Inc.

CS-100-1 v1.1



## Product Description

Order Number: **CS-100-1**

*Table 1. Contents and Storage Information*

Material	Amount	Storage
<b>Cleaning solution 1</b>	100 mL	RT

The detergent-based Cleaning Solution 1 (CS 1) is **ready-to-use**.

CS 1 is used for regular cleaning of the **heliX**<sup>cyto</sup> and **heliX**<sup>cyto</sup> chip. Please follow the use of CS 1 with RB 1 or another suitable buffer or H<sub>2</sub>O.

For research use only.

This product has a limited shelf life, please see expiry date on label.

## Contact

**Dynamic Biosensors GmbH**

Perchtinger Str. 8/10  
81379 Munich  
Germany

**Dynamic Biosensors, Inc.**

300 Trade Center, Suite 1400  
Woburn, MA 01801  
USA

**Order Information**   [order@dynamic-biosensors.com](mailto:order@dynamic-biosensors.com)

**Technical Support**   [support@dynamic-biosensors.com](mailto:support@dynamic-biosensors.com)

**[www.dynamic-biosensors.com](http://www.dynamic-biosensors.com)**

Instruments and chips are engineered and manufactured in Germany.

©2024 Dynamic Biosensors GmbH | Dynamic Biosensors, Inc. All rights reserved.