

Cleaning Solution 3

for maintenance of a heliXcyto

Dynamic Biosensors GmbH & Inc.

CS-50-3 v1.0



Product Description

Order Number: **CS-50-3**

Table 1. Contents and Storage Information

Material	Amount	Storage
Cleaning solution 3	50 mL	RT

Handle with personal protection (gloves, goggles) and avoid direct skin contact. Rinse thoroughly with plenty of water in case of skin contact.

Store in a dark, cool, and well aired location. Avoid contact with acids.

For research use only.

This product has a limited shelf life, please see expiry date on label.

Application Note

The Cleaning Solution 3 (CS 3) is ready-to-use for cleaning of a **helix^{cyto}** device on a **maintenance chip**.

DO NOT run over **helix^{cyto}** trapping chips!

Follow the use of CS 3 immediately with RB 1 or another suitable buffer or plenty of H₂O.

Contact

Dynamic Biosensors GmbH

Perchtinger Str. 8/10
81379 Munich
Germany

Dynamic Biosensors, Inc.

300 Trade Center, Suite 1400
Woburn, MA 01801
USA

Order Information order@dynamic-biosensors.com

Technical Support support@dynamic-biosensors.com

www.dynamic-biosensors.com

Instruments and chips are engineered and manufactured in Germany.

©2024 Dynamic Biosensors GmbH | Dynamic Biosensors, Inc. All rights reserved.