Article

# Programmable multispecific DNA-origami-based T-cell engagers

Received: 30 January 2023

Accepted: 26 June 2023

Published online: 17 August 2023

Check for updates

Klaus F. Wagenbauer<sup>1,2,7</sup>, Nhi Pham<sup>1,2,7</sup>, Adrian Gottschlich <sup>(3,4,7)</sup>, Benjamin Kick<sup>1,2,7</sup>, Viktorija Kozina<sup>1,2</sup>, Christopher Frank<sup>1,2</sup>, Daniela Trninic<sup>1,2</sup>, Pierre Stömmer <sup>(1,2)</sup>, Ruth Grünmeier <sup>(3)</sup>, Emanuele Carlini<sup>3</sup>, Christina Angeliki Tsiverioti<sup>3</sup>, Sebastian Kobold <sup>(3,5,6)</sup>, Jonas J. Funke<sup>1,2</sup> & Hendrik Dietz <sup>(1,2)</sup>

Multispecific antibodies have emerged as versatile therapeutic agents, and therefore, approaches to optimize and streamline their design and assembly are needed. Here we report on the modular and programmable assembly of IgG antibodies, F(ab) and scFv fragments on DNA origami nanocarriers. We screened 105 distinct quadruplet antibody variants in vitro for the ability to activate T cells in the presence of target cells. T-cell engagers were identified, which in vitro showed the specific and efficient T-cell-mediated lysis of five distinct target cell lines. We used these T-cell engagers to target and lyse tumour cells in vivo in a xenograft mouse tumour model. Our approach enables the rapid generation, screening and testing of bi- and multispecific antibodies to facilitate preclinical pharmaceutical development from in vitro discovery to in vivo proof of concept.

Programmable self-assembly with DNA origami enables fabricating discrete nanoscale objects with structurally well-defined two-dimensional and three-dimensional shapes from DNA molecules<sup>1–5</sup>, including nanoscale devices<sup>6–8</sup>, functional materials<sup>9,10</sup> and higher-order objects<sup>11,12</sup>. DNA origami objects are addressable and can be modified with various biomolecules in a site-specific fashion<sup>13,14</sup>. Previous studies have demonstrated the attachment of antibodies to DNA origami objects<sup>15,16</sup>, the binding of antibody-conjugated DNA origami objects to cell surfaces<sup>14,17,18</sup> and the modulation of T-cell function<sup>19,20</sup>. Recent developments such as the cost-efficient mass production of DNA origami raw materials and stabilization approaches for in vivo application<sup>21-23</sup> may enable the clinical translation of diverse therapeutical concepts such as DNA origami biomedical nanorobots<sup>10,24</sup>.

In parallel to the advances in DNA nanotechnology, cancer immunotherapies have contributed to a paradigm shift in oncological treatment landscapes<sup>25</sup>. In particular, T-cell-centred immunotherapies (for example, immune-checkpoint-inhibiting antibodies) are now established in clinical practice in various cancer entities<sup>26</sup>. In addition, in B-cell-derived haematological malignancies (such as acute lymphoblastic leukaemia (ALL) or B-cell lymphomas), T-cell-engaging antibodies have led to clinical responses even in otherwise treatment-refractory patients<sup>27</sup>. The Federal Drug Agency (FDA) and European Medicine Agency (EMA) have granted approval for blinatumomab, a CD19-CD3-bispecific T-cell engager (BiTe), prolonging the overall survival of patients.

In B-cell malignancies, using B-cell-associated antigens such as CD19 and CD20 has proven feasible, efficacious and manageable from a safety perspective, partly owing to established clinical treatments to manage induced B-cell aplasia<sup>28</sup>. However, this is unique to B-cell-targeting agents and cannot be expected in other diseases. In solid cancers, tumour-associated antigens are often co-expressed on vital epithelial tissues, creating the risk for severe on-target off-tumour

<sup>1</sup>Department of Biosciences, School of Natural Sciences, Technical University of Munich, Garching, Germany. <sup>2</sup>Munich Institute of Biomedical Engineering, Technical University of Munich, Garching, Germany. <sup>3</sup>Division of Clinical Pharmacology, University Hospital, LMU Munich, Member of the German Center for Lung Research (DZL), Munich, Germany. <sup>4</sup>Department of Medicine III, University Hospital, LMU Munich, Munich, Germany. <sup>5</sup>German Cancer Consortium (DKTK), Partner Site Munich, Munich, Germany. <sup>6</sup>Einheit für Klinische Pharmakologie (EKLiP), Helmholtz Munich, Research Center for Environmental Health (HMGU), Neuherberg, Germany. <sup>7</sup>These authors contributed equally: Klaus F. Wagenbauer, Nhi Pham, Adrian Gottschlich, Benjamin Kick. e-mail: sebastian.kobold@med.uni-muenchen.de; jonasfunke@mytum.de; dietz@tum.de toxicities<sup>29</sup>. Increasing cell-type specificity, for example, by the simultaneous targeting of multiple-tumour-associated antigens has the potential to minimize the risk for severe on-target off-tumour toxicities<sup>30</sup>. In addition, targeting more than one antigen on a target cell may prove beneficial in preventing antigen-negative relapse, as sequential or simultaneous multiple targeting will enhance therapeutic pressure and counter the development of negative variants.

Thus, multispecific molecules targeting cancer vulnerabilities are needed to leverage immune cell potential in oncology. A large number of drug candidates are currently in preclinical and clinical development, with the focus shifting from bispecific antibodies and BiTe formats (four on the market and more than 100 in clinical development) towards formats with increased specificities or enhanced pharmacokinetic properties (eight candidates in clinical development)<sup>31-33</sup>. Various approaches have been developed to produce multispecific antibodies, most of which rely on fusing engineered antibody domains<sup>34</sup>. Although these approaches support controlling the degree of valency, the spatial geometric arrangements of the domains are restricted by the structural constraints of the protein scaffold used.

Here we used programmable self-assembly with DNA origami to create a synthetic antibody carrier platform called programmable T-cell engager (PTE). PTEs offer desirable properties for T-cell engagement, including the capability to modularly position antibodies (IgG, F(ab) or scFv) with control over valency, orientation and spatial arrangement. We provide the proof of concept of specific T-cell engaging in vitro and validate the PTE functionality in leukaemia models in vivo.

#### Assembly and screening of IgG-based PTEs

The ability to place IgG antibodies in a user-defined fashion on a DNA origami carrier is the prerequisite for building more complex multivalent configurations. We tested and optimized methods to meet these requirements in auxiliary experiments (Supplementary Notes 1 and 2 and Supplementary Figs. 1-3), including demonstrating multivalent cell binding and testing for cell internalization (Supplementary Note 3 and Supplementary Figs. 4 and 5). Building on these optimized methods, we created a tetravalent antibody carrier featuring four distinct antibody attachment sites (Fig. 1a). For each attachment site on the DNA origami chassis, we created a library of DNA-tagged antibodies carrying a sequence-complementary single-stranded DNA tag. To induce the activation of effector T cells, we chose anti-CD3, anti-CD28 and anti-CD137 antibodies (Fig. 1a, left). To mediate binding to the target cells, we chose antibodies against the known antigens CD19, CLL-1, CD22 and CD123. We prepared antibody-DNA conjugates and then assembled 105 unique antibody combinations. We validated the assembly of these PTEs via gel electrophoretic mobility analysis (Fig. 1b) and quantified the yield of fully assembled tetravalent combinations for each variant. The assembly yield varied between from 85% (variant 2× anti-CD123 2× anti-CD3) to 97% (variant 2× anti-CD19 2× anti-CD3).

To study PTE-induced T-cell activation, we used a nuclear factor of activated T cells (NFAT)-luciferase reporter assay<sup>35</sup>. We co-cultured CD19<sup>+</sup>, CD22<sup>+</sup>, CD123<sup>-</sup> and CLL-1<sup>-</sup> NALM-6 ALL cells with NFAT-luciferase-transduced Jurkat cells in the presence or absence of PTEs (Fig. 1c,d). Bivalent aCD3 IgG antibodies can crosslink T-cell receptors causing T-cell activation in the absence of target cells. We, therefore, subtracted the background signal generated by the PTEs that carried only T-cell antibodies (Supplementary Fig. 6). Here aCD19× aCD3 constructs lead to the strongest T-cell activation, whereas aCD22 × aCD3 constructs induced weak T-cell activation (Fig. 1c). These results are in line with the differing densities of CD19 and CD22 molecules on NALM-6 cells<sup>36</sup>. Also, aCLL-1, aCD123 × aCD3 constructs did not cause target-cell-induced T-cell activation, proving the specificity of the PTEs. The inclusion of additional T-cell-activating antibodies (aCD28, aCD137) resulted in the increased activation of T cells compared with the single aCD19 × aCD3 variant, with the aCD19 × aCD3-aCD28 PTEs showing the strongest T-cell activation (Fig. 1c). Again, these results are compliant with the literature<sup>37</sup>.

Next, we analysed the dual-tumour-targeting constructs for their ability to activate T cells (Fig. 1d). We observed increased activation signals for two-target variants with a single activation of T cells via CD3 compared with single-target variants (Fig. 1d). The addition of a co-stimulatory domain to the dual-targeting variants further enhanced the signal. We observed the strongest activation for variants consisting of aCD3, aCD28 and two aCD19 antibodies or one aCD19 antibody and one aCD22 antibody (Fig. 1e). The aCLL-1 control constructs did not induce the activation of T cells.

#### In vitro characterization of F(ab)-based PTEs

Protein-based bispecific antibodies can mediate the potent T-cell-directed lysis of tumour cells<sup>38</sup>. Our approach should allow us to assemble bispecific variants with similar capabilities. However, the size of the antibody and the arrangement of the paratopes are essential factors that can influence the efficacy of target cell killing<sup>39</sup> (Supplementary Fig. 7). In this context, full-sized IgGs present limitations, such as Fc-domain-mediated binding to immune cells, which can cause undesired cell interactions<sup>40</sup>, and crosslinking receptors through their two paratopes, leading to non-specific target cell lysis (Supplementary Fig. 8). With these aspects in mind, we designed a smaller antibody carrier chassis  $(20.0 \times 15.0 \times 7.5 \text{ nm}^3)$ that can display an anti-CD3 F(ab) fragment for T-cell binding and up to four F(ab) fragments to recognize the target cell antigens (Fig. 2a,b and Supplementary Fig. 9). We placed the F(ab) fragments at the corners of the DNA origami chassis to realize a paratope-to-paratope distance of approximately 7.5 nm, which is comparable with reported cell-cell distances in immunological synapses<sup>41</sup>. The antibody attachment concepts established with the larger chassis were directly transferable to the small chassis and enabled fabricating multispecific variants with high yields (>98%), as seen by agarose gel electrophoresis and transmission electron microscopy (TEM) (Fig. 2b,c and Supplementary Fig. 10). To orient the F(ab) fragments on the DNA origami chassis, we site-specifically coupled the adapter DNA strands to one of the thiol groups of its disulfide bond, located opposite the paratope (Supplementary Fig. 11).

Using the smaller chassis, we prepared five bispecific PTE variants targeting four known human-tumour-associated antigens, namely, huCD19 (PTE-2×19-3), huCD33 (PTE-2×33-3), huEpCAM (PTE-2×EpCAM-3) and huPSMA (PTE-2×PSMA-3) or murine CD19 (mPTE-2×mCD19-3) (Fig. 2d-h). First, we tested if PTE-2×19-3 constructs induced the lysis of CD19<sup>+</sup> NALM-6 cells in co-cultures with peripheral blood mononuclear cells (PBMC). Using flow cytometry, we observed the concentration-dependent lysis of the NALM-6 target cells after 24 h of co-culture. The lysis efficiency was higher than 90% at concentrations above 500 pM with a calculated potency of approximately 100 pM (EC<sub>50</sub> value; Supplementary Fig. 12). In contrast, the monospecific controls (small chassis with one ahuCD3 (PTE-3) or small chassis with ahuCD19 only (PTE-19)) induced minimal target cell lysis, corroborating the specificity of the designed constructs (Fig. 2d). In parallel, we also monitored the activation of CD8<sup>+</sup>T cells by analysing the CD69 expression (Supplementary Fig. 13b). Again, only the bispecific construct activated the T cells, whereas PTE-3 only showed low background activation at high concentrations. The lysis efficacy and T-cell activation efficiency of our small chassis variant at 1 nM concentration and above were comparable with an huCD19 × huCD3 BiTe (blinatumomab-biosimilar (Blina-BS), InvivoGen). Blina-BS showed a superior lysis efficacy at concentrations below 0.1 nM relative to our compounds. We attribute the potency differences primarily to different affinities of the antibody paratopes used in our construct relative to those in the BiTe molecule. To test the transferability of the concept to other antigens, we tested the PTE variants against huCD33, huEpCAM, huPSMA and muCD19 in cytotoxic T-cell assays, with tumour cells expressing the respective antigens. To this end, PTE-2×33-3 were tested in co-cultures with acute myeloid leukaemia cell line CD33<sup>+</sup> MOLM-13, whereas PTE-2×PSMA-3



**Fig. 1** | **Production and functional screening of 105 unique antibody combinations on a DNA chassis. a**, Schematic of a multispecific antibody chassis variant library created from a set of antibody–DNA conjugates. The symbols indicate the antibody, and the colour indicates the engaged cell type. Antibodies are covalently tagged with DNA handles with the sequences A, B, C or D, depending on the library, and the sequences are complementary to DNA handles on the chassis (centre). The chassis carries four DNA handles. Antibody chassis variants are produced by mixing the respective antibodies from the libraries with the DNA chassis. Variants are named by their antibody combination (the centre bottom shows an example combination). Two reference-free class averages calculated from single-particle TEM micrographs. Scale bar, 20 nm. The top average shows the platform without antibodies and the bottom average image shows the platform with four IgG antibodies, indicated by the orange and

blue arrow heads pointing to the blurred additional signal in the average image. b, Montage of laser-scanned images of agarose gels on which 105 variants were electrophoresed that were incubated with different antibody combinations (as indicated by the symbols). The first and last lanes show a reference 5 MDa DNA origami object. c,d, T-cell activation was measured by using NFAT-luciferase Jurkat cell line in co-cultures with human ALL cell line NALM-6 in the presence of the indicated combinations (c,d). Relative T-cell activation (normalized to variants without target cell antibodies) of different variants for 100 pM and 1,000 pM DNA chassis concentrations. The icons in orange and blue indicate the respective antibodies used in the combination. e, Relative T-cell activation of the variants sorted for maximum activation according to the values in c and d at 100 pM.

or PTE-2×EpCAM-3 were analysed in co-cultures with PSMA<sup>+</sup> LNCaP prostate cancer or EpCAM<sup>+</sup> MCF-7 breast cancer cells, respectively (Fig. 2e–g). All the PTE variants induced the dose-dependent lysis of target cells with high lysis efficacy in the low nanomolar range. To evaluate the transferability of our approach to murine surrogate molecules, mPTE-3×mCD19-3 were tested in the co-cultures of murine splenocytes with mCD19<sup>+</sup> A20 cells (Fig. 2h). Using live-cell imaging, we directly observed the specific T-cell-mediated lysis of A20 cells for the bispecific molecule (Fig. 2i).

Finally, we used the smaller F(ab)-based chassis to evaluate our initial findings with the larger IgG-based chassis. To study how the target-cell-binding valency impacts the lysis efficacy of our smaller chassis, we assembled and tested the PTE variants with up to four ahuCD19 F(ab) fragments and one ahuCD3 F(ab) fragment in cyto-toxic T-cell assays (Supplementary Fig. 14). The potency increased substantially from one to two aCD19 F(ab) fragments did not improve the potency.



**Fig. 2** | **Programmable T-cell-mediated killing of target cells. a**, Schematic of a multispecific brick-shaped antibody carrier (chassis) with dimensions of 10.0 × 15.0 × 7.5 nm<sup>3</sup>. The grey cylinders represent DNA double helices, and the F(ab) fragments are coloured in orange (anti-CD3) and blue (representing a F(ab) fragment for antigens located on the target cells). **b**, Negatively stained TEM image of the small chassis (left) and negatively stained TEM image of the small chassis (left) and negatively stained TEM image of the small chassis with 1× anti-CD3 and 2× anti-CD19 F(ab) fragments (right). Scale bar, 25 nm. The arrows in blue and yellow highlight the attached F(ab) fragments as an example. **c**, Laser-scanned image of an agarose gel on which different samples were electrophoresed. The samples were prepared with different F(ab) fragment combinations (as indicated by the numbers). P, pocket; icons highlight the different antibody chassis variants; 0-0, reference for the migration of platform only. **d**-**h**, Cytotoxic T-cell-mediated target cell lysis assays. Fraction

of dead target cells after 24 h as a function of PTE concentration in the assay (Supplementary Information). Effector (PBMC) and target cell ratio was chosen as 5:1. Red dots, multispecific T-cell-engaging variant with anti-CD3 and at least two target-specific F(ab) fragments for the respective cell line. Black dots, monospecific controls. Solid lines, Hill fit to the data. Dashed line, PBMC and target cells without PTE after 24 h. Grey squares, blinatumomab-biosimiliar (Blina-BS). The error bars to the data are standard deviations to the mean of three biological replicates. **i**, Live-cell fluorescence microscopy over 24 h of a mixture containing A20 cells (stained with CellTrace CSFE; blue) and splenocytes in a 1:5 ratio, and a variant (1× mu anti-CD3-4× mu anti-CD19) carrying a fluorescent tag (cyanine-5). A live-dead stain was used for the visualization of dead cells (SYTOX Orange; cyan). Scale bar, 2.5 µm.

#### **Biodistribution of PTEs**

To ensure PTE stability under physiological low-ionic-strength conditions and achieve protection against nuclease degradation, PTEs were coated with PEG-oligolysine<sup>22</sup>. We validated the successful assembly of PTEs in combination with PEG-oligolysine coating via mass photometry (Supplementary Fig. 15). These PTEs showed storage stability at 4 °C over at least three months with no sign of degradation or loss of functionality in T-cell-mediated lysis assays (Supplementary Fig. 16). Incubation in fetal bovine serum (FBS) (100% serum concentration) for four hours followed by a 24 h T-cell-mediated lysis assay (50% serum concentration) of NALM-6 cells showed comparable results to samples incubated with heat-inactivated FBS (Supplementary Fig. 17). In contrast, PTE samples without a coating showed a substantial decrease in lysis efficiency. We also analysed the composition of the solvent regarding impurities such as endotoxins that may cause immunological side effects in vivo. All the administered samples were below 36 EU ml<sup>-1</sup> for a daily dose of 100 µl, satisfying the standard residual concentration of endotoxins for in vivo applications (Supplementary Methods)42.

We used in vivo time-resolved fluorescence imaging to measure the distribution of PTEs in tumour-free immunodeficient NOD-SCID IL-2Ry-null (NSG) mice. Mice received one intravenous (i.v.) injection (15 pmol) of either phosphate-buffered saline (PBS), DNA chassis or PTE-2×19-3, with fluorescent Cy7 dyes incorporated into the DNA chassis or a mixture of 20 Cy7-modified DNA strands only (here termed Cy7-strand) (Fig. 3a). After injection of the samples, we measured Cy7 fluorescence in the anaesthetized state every 30 s for a duration of 2.5 h (Fig. 3b). Mice treated with the Cy7-dye-containing samples (DNA chassis, PTE-2×19-3 and Cy7-strand) had an increased fluorescence radiance signal compared with the vehicle control.

Five minutes after i.v. injection, the Cy7 signal was delocalized across the mice with a slight accumulation at the liver for all the samples. The fluorescence signal then gradually localized to the liver and to the bladder regions. Accumulation in the bladder region was observed more quickly in the mice that received Cy7-strands or DNA chassis controls compared with PTE-2×19-3-treated mice. We analysed the Cy7 intensity of the liver, bladder and periphery versus time for each variant (Fig. 3c). Overall, traces indicate that the Cy7-strands and DNA chassis controls accumulate in the bladder and the liver regions, whereas PTE-2×19-3 mainly localized in the liver region. Mice given the Cy7 control treatment cleared out faster at the periphery than those administered PTE-2×19-3. We analysed the distribution of the samples



**a**, Schematic of the experimental procedure. NSG mice were administered the highlighted samples carrying a Cy7 dye via a tail-vein injection. Distribution of the administered sample was then continuously measured using an in vivo imaging platform (IVIS). Colour scheme, bispecific PTE with Cy7 (midnight; PTE-2×19-3), DNA oligo with Cy7 (mocha; Cy7 DNA), DNA origami chassis with Cy7 (dark red; Cy7 chassis) and vehicle (grey; PBS). **b**, IVIS fluorescent images of NSG mice over

time. Left to right, samples as indicated in **a**. **c**, Fluorescence intensity of the bispecific (bluish), Cy7 DNA (brownish) and Cy7 chassis (reddish) samples as a function of time across different organs/regions. Data are mean  $\pm$  standard error of mean (s.e.m.) from three independent measurements of n = 3 mice. **d**, **e**, IVIS fluorescence images of different organs and excreta extracted from n = 3 NSG mice. **f**, **g**, Fluorescence intensity in different organs and excreta. The background fluorescence of each organ was subtracted. Data are mean  $\pm$  s.e.m. from n = 3 mice.

in the harvested organs at the end of the experiment (Fig. 3d–g). In accordance with the in vivo imaging results, the strongest Cy7 signal was detected in the colon, liver and kidneys (Fig. 3e,g). However, the Cy7 signal was also measurable in the brain, heart, lungs, spleen and bone marrow (Fig. 3d,f). These results are consistent with previous studies<sup>43,44</sup> and indicate that DNA-origami-based PTEs readily distribute in living animals and are primarily eliminated through the biliary and renal excretory systems.

#### In vivo characterization of PTEs

To target T cells against tumour cells in living organisms, PTEs must specifically bind to and remain bound to the target cells until a T cell is recruited. To investigate these processes, we intravenously injected NSG mice first with NALM-6-GFP-luciferase (luc) target cells, then with PBMC effector cells and finally with a single dose of vehicle (PBS), Cy5-modified PTE variant (PTE-2×19-3, PTE-3 and PTE-2×19) or Blina-BS (Fig. 4a). After 4 h, the mice were sacrificed and cells from the



**Fig. 4** | **Investigation of functionality of DNA-origami-based T-cell engagers in vivo. a**, Treatment schedule of binding and activation trials. Experimental procedure to measure the in vivo binding and activation behaviour on tumour (NALM-6) and T cells. Colour scheme: midnight blue, PTE-2×19-3; teal, PTE-3; mocha, PTE-2×19; red, chassis; violet, Blina-BS. Doses were between 10 and 15 pmol. NALM-6-GFP-luc, PBMC and samples were intravenously administered at the indicated time points. Three mice were injected per group. After 4 h, the mice were sacrificed and sample distribution in the bone marrow was analysed using flow cytometry. **b,c**, Mean fluorescence intensity (MFI) measurements (left) and representative histograms (right) of flow cytometric measurements of Cy5 fluorescence (**b**) or CD19-BV785 fluorescence (**c**) on NALM-6-GFP-luc cells (indicating Cy5-labelled PTE bound to the tumour cells). **d**, **e**, Measurement of CD69 fluorescence on either CD4<sup>+</sup> (**d**) or CD8<sup>+</sup> (**e**) T cells. In **b**–**e**, data are mean  $\pm$  s.e.m. from n = 3 mice. **f**, Schematic of the long-term treatment with PTEs in vivo. **g**, Luminescence images of mice injected with luciferin to visualize NALM-6-GFP-luc tumour cells. The white cross indicates the mouse that was censored due to non-tumour-related toxicities. **h**, Quantification of bioluminescence measurements. For **g** and **h**, n = 5 mice per group; for **h**, mean  $\pm$  s.e.m. is shown. For all the panels, statistical significance was calculated using ordinary one-way or two-way analysis of variance with Tukey's multiple comparisons correction.

bone marrow were analysed using flow cytometry (Fig. 4b–e).  $GFP^+$ NALM-6 target cells showed an increased Cy5 fluorescence signal only in mice receiving PTE variants that carry ahuCD19 F(ab) fragments (PTE-2×19-3 and PTE-2×19) (Fig. 4b). In addition, we probed the CD19 antigen occupancy by staining the cells with an anti-CD19 antibody before flow cytometric analysis. NALM-6 cells from the vehicle and PTE-3 samples had an increased CD19 signal, indicating the accessibility of the CD19 antigen (Fig. 4c). In contrast, the signal of the anti-CD19 antibody was reduced in mice treated with PTE-2×19-3, PTE-2×19 or Blina-BS (Fig. 4c), indicating that CD19 epitopes are occupied by either PTEs (76% occupied) or Blina-BS (99% occupied), respectively. Overall, these results suggest that PTEs specifically bind to the target cells in vivo and remain bound for several hours (>4 h). As demonstrated in vitro, by crosslinking effector cells and target cells, T-cell engagers recruit and activate the T cells. To validate this mechanism in vivo, we quantified T-cell activation by measuring the expression of CD69 on the transferred T cells as a marker for early T-cell activation (Fig. 4d,e). Both Blina-BS and PTE-2×19-3 induced a significant increase in CD69 expression on CD4<sup>+</sup> or CD8<sup>+</sup> T cells compared with the control constructs (Fig. 4d,e). Our results, therefore, indicate that PTEs bind NALM-6 target cells and activate the T cells in vivo.

Next, we studied whether PTEs can also control the tumour outgrowth of NALM-6 target cells in vivo. To this end, we first determined the optimal PTE dose. We administered 20 million NALM-6-GFP-luc cells intravenously and then treated the mice with different concentrations of PTE-2×19-3 (300, 100 and 30 pmol) or control samples (PTE-3, vehicle and Blina-BS) of equivalent doses (Supplementary Fig. 19a). PBMC were intravenously injected on the days of treatment. Bioluminescence (BLI) measurements after six days revealed a significant reduced tumour burden in mice that received PTE-2×19-3 compared with PTE-3- or vehicle-treated mice (Supplementary Fig. 19b,c). Flow cytometric analysis of the bone marrow confirmed the reduced tumour burden for PTE-2×19-3-treated mice and revealed a dose-dependent tumour reduction (Supplementary Fig. 19c). To determine the optimal anti-CD3 clone, we compared the T-cell activation of PTE variants with OKT3 and UCHT1 CD3 binders (Supplementary Fig. 19d). Again, NALM-6-GFP-luc cells were intravenously administered and we injected PBMC and different concentrations of a UCHT1-based 2×19-3 PTE (300 and 30 pmol) or OKT3-based 2×19-3 PTE (30 pmol) or control samples (vehicle and Blina-BS). Flow cytometric analysis of the bone marrow revealed the activation of T cells in all mice that received PTE-2×19-3, with minimal differences between the different CD3 binders (Supplementary Fig. 19e, f). However, mice treated with UCHT1-based PTEs exhibited dose-dependent T-cell depletion compared with mice treated with OKT3-based PTEs (Supplementary Fig. 20).

We then investigated whether OKT3-based PTEs can control tumour outgrowth for an extended duration. We again intravenously injected NALM-6-GFP-luc tumour cells into NSG mice. After three days, mice repeatedly received PBMC and PTE-2×19-3 or PTE-3 and Blina-BS as negative or positive controls, respectively. Tumour growth was subsequently monitored with BLI for a total of 21 days (Fig. 4f). PTE-2×19-3 or Blina-BS treatment led to reduced tumour growth in the treated mice, compared with mice that received PTE-3-negative control samples (Fig. 4g,h). In summary, our studies demonstrate the functionality of PTEs built from DNA origami chassis in vivo.

#### Conclusion

DNA origami allows the precise spatial arrangement of biomolecules, and when combined with antibodies, this approach represents an exciting avenue for developing biomedical nanodevices. Here we created a programmable DNA origami chassis that allows the positioning of multiple antibodies and demonstrated the functionality both in vitro and in vivo. By optimizing the attachment strategy as well as the purification and stabilization methods, we demonstrated the ability to assemble and screen a large number of variants in parallel. For example, we assembled 105 multispecific antibody variants to identify and rank antibody combinations for T-cell activation. Such PTEs can target multiple copies of the same antigen and also sets of distinct antigens, allowing for programming more advanced cells or antigen recognition.

The ability to display custom sets of antibodies on DNA origami chassis allowed us to design PTEs that display antibodies against tumour cells and T cells, thereby directing T cells to specifically lyse tumour cells. Although our PTEs and previously developed protein-based aCD19-aCD3 antibodies have comparable efficacies, the PTEs at present require higher concentrations to achieve the same biological effect<sup>45</sup>. The potency differences may be attributed to, for example, different binding modes of the antibody clones used, not-vet-optimal placement of F(ab) fragments on the DNA origami chassis with respect to redirecting T-cell lysis and cellular internalization of the DNA origami PTEs competing with target cell recruitment and lysis<sup>46,47</sup>. The approach may be extended in the future to simultaneously engage several signalling pathways for improved T-cell activation in one molecule. For example, we have also tested combining ahuCD3 with ahuCD28 antibodies (Fig. 1 and Supplementary Fig. 21). The addition of ahuCD28 caused a substantially increased T-cell activation signal for IgG-based PTEs, but only slightly enhanced the potency in our in vitro lysis assays with F(ab)-based PTEs. CD28 co-stimulation may be more beneficial for other tumour models, where the microenvironment and the exhaustion of T cells are important factors for tumour lysis efficacy.

PTEs functioned in vivo and distributed well within animals<sup>22,48</sup>. We showed that our PTEs specifically bind to target cells and recruit T cells in vivo, thus demonstrating the mechanism of action within living animals. In efficacy experiments, our PTEs control tumour outgrowth, demonstrating the applicability of DNA-origami-based therapeutics for cancer treatment.

Given the DNA origami technologies' modularity, adjustability and high degree of addressability, we expect that a wide range of complex and even logic-gated chassis<sup>10</sup> for immunotherapy can be engineered. Such chassis have the potential to help overcome important challenges currently faced in the field, including discriminating healthy cells from tumour cells based on surface markers by detecting patterns of multiple antigens rather than on single targets. Additionally, our DNA origami chassis enables the simultaneous engagement of several hallmark signalling pathways for improved T-cell activation in one molecule, including checkpoint inhibition but also targeted co-stimulation. We believe our results will enable the clinical application of DNA nanotechnologies and highlight the potential of DNA-origami-based biomolecular engineering strategies for medical applications.

#### **Online content**

Any methods, additional references, Nature Portfolio reporting summaries, source data, extended data, supplementary information, acknowledgements, peer review information; details of author contributions and competing interests; and statements of data and code availability are available at https://doi.org/10.1038/s41565-023-01471-7.

#### References

- 1. Rothemund, P. W. K. Folding DNA to create nanoscale shapes and patterns. *Nature* **440**, 297–302 (2006).
- 2. Douglas, S. M. et al. Self-assembly of DNA into nanoscale three-dimensional shapes. *Nature* **459**, 414–418 (2009).
- 3. Han, D. et al. DNA origami with complex curvatures in three-dimensional space. *Science* **332**, 342–346 (2011).
- 4. Benson, E. et al. DNA rendering of polyhedral meshes at the nanoscale. *Nature* **523**, 441–444 (2015).
- 5. Dietz, H., Douglas, S. M. & Shih, W. M. Folding DNA into twisted and curved nanoscale shapes. *Science* **325**, 725–730 (2009).
- 6. Funke, J. J. & Dietz, H. Placing molecules with Bohr radius resolution using DNA origami. *Nat. Nanotechnol.* **11**, 47–52 (2016).
- Zhang, Z., Yang, Y., Pincet, F., Llaguno, M. C. & Lin, C. Placing and shaping liposomes with reconfigurable DNA nanocages. *Nat. Chem.* 9, 653–659 (2017).
- 8. Knudsen, J. B. et al. Routing of individual polymers in designed patterns. *Nat. Nanotechnol.* **10**, 892–898 (2015).
- Auvinen, H. et al. Protein coating of DNA nanostructures for enhanced stability and immunocompatibility. *Adv. Healthc. Mater.* 6, 1700692 (2017).

- Douglas, S. M., Bachelet, I. & Church, G. M. A logic-gated nanorobot for targeted transport of molecular payloads. *Science* 335, 831–834 (2012).
- Tikhomirov, G., Petersen, P. & Qian, L. Fractal assembly of micrometre-scale DNA origami arrays with arbitrary patterns. *Nature* 552, 67–71 (2017).
- Wagenbauer, K. F., Sigl, C. & Dietz, H. Gigadalton-scale shape-programmable DNA assemblies. *Nature* 552, 78–83 (2017).
   Original Mathematical Mathematical Mathematical Activity of the statement of the
- Seitz, I., Ijäs, H., Linko, V. & Kostiainen, M. A. Optically responsive protein coating of DNA origami for triggered antigen targeting. ACS Appl. Mater. Interfaces 14, 38515–38524 (2022).
- Cremers, G. A. et al. Determinants of ligand-functionalized DNA nanostructure-cell interactions. J. Am. Chem. Soc. 143, 10131–10142 (2021).
- Shaw, A. et al. Binding to nanopatterned antigens is dominated by the spatial tolerance of antibodies. *Nat. Nanotechnol.* 14, 184–190 (2019).
- 16. Zhang, P. et al. Capturing transient antibody conformations with DNA origami epitopes. *Nat. Commun.* **11**, 3114 (2020).
- Rosier, B. J. et al. Incorporation of native antibodies and Fc-fusion proteins on DNA nanostructures via a modular conjugation strategy. *Chem. Commun.* 53, 7393–7396 (2017).
- 18. Akbari, E. et al. Engineering cell surface function with DNA origami. *Adv. Mater.* **29**, 1703632 (2017).
- Sun, Y. et al. DNA origami based artificial antigen-presenting cells for adoptive T cell therapy. Sci. Adv. 8, eadd1106 (2022).
- 20. Fang, T. et al. Spatial regulation of T-cell signaling by programmed death-ligand 1 on wireframe DNA origami flat sheets. *ACS Nano* **15**, 3441–3452 (2021).
- 21. Praetorius, F. et al. Biotechnological mass production of DNA origami. *Nature* **552**, 84–87 (2017).
- 22. Ponnuswamy, N. et al. Oligolysine-based coating protects DNA nanostructures from low-salt denaturation and nuclease degradation. *Nat. Commun.* **8**, 15654 (2017).
- Gerling, T., Kube, M., Kick, B. & Dietz, H. Sequence-programmable covalent bonding of designed DNA assemblies. *Sci. Adv.* 4, eaau1157 (2018).
- 24. Li, S. et al. A DNA nanorobot functions as a cancer therapeutic in response to a molecular trigger in vivo. *Nat. Biotechnol.* **36**, 258–264 (2018).
- 25. Kruger, S. et al. Advances in cancer immunotherapy 2019 latest trends. *J. Exp. Clin. Cancer Res.* **38**, 268 (2019).
- Yu, J. X., Hubbard-Lucey, V. M. & Tang, J. Immuno-oncology drug development goes global. *Nat. Rev. Drug Discov.* 18, 899–901 (2019).
- 27. Wang, S. et al. The state of the art of bispecific antibodies for treating human malignancies. *EMBO Mol. Med.* **13**, e14291 (2021).
- Topp, M. & Feuchtinger, T. Management of hypogammaglobulinaemia and B-cell aplasia. in *The EBMT/EHA* CAR-T Cell Handbook 147–149 (Springer, 2022).
- Tokarew, N., Ogonek, J., Endres, S., von Bergwelt-Baildon, M. & Kobold, S. Teaching an old dog new tricks: next-generation CAR T cells. *Br. J. Cancer* **120**, 26–37 (2019).
- Choe, J. H. et al. SynNotch-CAR T cells overcome challenges of specificity, heterogeneity, and persistence in treating glioblastoma. Sci. Transl. Med. 13, eabe7378 (2021).
- Banaszek, A. et al. On-target restoration of a split T cell-engaging antibody for precision immunotherapy. *Nat. Commun.* **10**, 5387 (2019).
- 32. Geis, M. et al. Combinatorial targeting of multiple myeloma by complementing T cell engaging antibody fragments. *Commun. Biol.* **4**, 44 (2021).
- 33. Ma, J. et al. Bispecific antibodies: from research to clinical application. *Front. Immunol.* **12**, 626616 (2021).

- Elshiaty, M., Schindler, H. & Christopoulos, P. Principles and current clinical landscape of multispecific antibodies against cancer. *Int. J. Mol. Sci.* 22, 5632 (2021).
- Stecha, P. et al. Improved T Cell Activation Bioassay for Development of Bispecific Antibodies and Engineered T Cell Immunotherapies (Promega, 2016).
- Haso, W. et al. Anti-CD22–chimeric antigen receptors targeting B-cell precursor acute lymphoblastic leukemia. *Blood* 121, 1165–1174 (2013).
- Trickett, A. & Kwan, Y. L. T cell stimulation and expansion using anti-CD3/CD28 beads. J. Immunol. Methods 275, 251–255 (2003).
- Baeuerle, P. A. & Reinhardt, C. Bispecific T-cell engaging antibodies for cancer therapy. *Cancer Res.* 69, 4941–4944 (2009).
- Bluemel, C. et al. Epitope distance to the target cell membrane and antigen size determine the potency of T cell-mediated lysis by BiTE antibodies specific for a large melanoma surface antigen. *Cancer Immunol. Immunother.* 59, 1197–1209 (2010).
- 40. Junker, F., Gordon, J. & Qureshi, O. Fc gamma receptors and their role in antigen uptake, presentation, and T cell activation. *Front. Immunol.* **11**, 1393 (2020).
- 41. McCann, F. E. et al. The size of the synaptic cleft and distinct distributions of filamentous actin, ezrin, CD43, and CD45 at activating and inhibitory human NK cell immune synapses 12. *J. Immunol.* **170**, 2862–2870 (2003).
- 42. Malyala, P. & Singh, M. Endotoxin limits in formulations for preclinical research. J. Pharm. Sci. **97**, 2041–2044 (2008).
- 43. Wamhoff, E.-C. et al. Evaluation of nonmodified wireframe DNA origami for acute toxicity and biodistribution in mice. *ACS Appl. Bio Mater.* **6**, 1960–1969 (2023).
- 44. Mou, Q. et al. DNA trojan horses: self-assembled floxuridine-containing DNA polyhedra for cancer therapy. *Angew. Chem. Int. Ed.* **56**, 12528–12532 (2017).
- Löffler, A. et al. A recombinant bispecific single-chain antibody, CD19×CD3, induces rapid and high lymphoma-directed cytotoxicity by unstimulated T lymphocytes. *Blood* 95, 2098–2103 (2000).
- Wang, P. et al. Visualization of the cellular uptake and trafficking of DNA origami nanostructures in cancer cells. *J. Am. Chem. Soc.* 140, 2478–2484 (2018).
- Bastings, M. M. C. et al. Modulation of the cellular uptake of DNA origami through control over mass and shape. *Nano Lett.* 18, 3557–3564 (2018).
- 48. Balakrishnan, D., Wilkens, G. D. & Heddle, J. G. Delivering DNA origami to cells. *Nanomedicine* **14**, 911–925 (2019).

**Publisher's note** Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

**Open Access** This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons license, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons license and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this license, visit http://creativecommons. org/licenses/by/4.0/.

 $\circledast$  The Author(s), under exclusive licence to Springer Nature Limited 2023

#### Article

#### Methods

All the experiments and experimental conditions described throughout this study comply with the ethical regulations set forth by the institutional review board of the medical faculty of the Ludwig-Maximilians-Universität and the Regierung von Oberbayern (approval of animal experiments).

#### Chemicals, antibodies and cell lines

Unless otherwise mentioned, chemicals used within this work were purchased from Sigma-Aldrich and all the IgG antibodies, RPMI, PBS, FBS and penicillin–streptomycin were purchased from Thermo Fisher.

The human cell lines Jurkat (T-cell leukaemia, DSMZ, no. ACC-282), NALM-6 (B-cell precursor leukaemia, DSMZ, no. ACC-128), MCF-7 (breast adenocarcinoma, DSMZ, no. ACC-115) and Molm-13 (acute myeloid leukaemia, DSMZ, no. ACC 554) were obtained from DMSZ. Human cell line LNCaP (metastatic lesion of prostate adenocarcinoma, CLS, 300265) was obtained from CLS. All the cell lines were stored in liquid nitrogen. Successful cell-line authentication was done via polymerase chain reaction for Jurkat, NALM-6, LNCap and MCF-7.

Jurkat, NALM-6 and Molm-13 cells were grown in RPMI 1640 medium, supplemented with 10% FBS, penicillin (200 U ml<sup>-1</sup>), streptomycin (200  $\mu$ g ml<sup>-1</sup>) and with or without additional 20 mML-glutamine for Jurkat and NALM-6, respectively. The MCF-7 cells were grown in high-glucose (25 mM D-glucose) Dulbecco's modified Eagle's medium medium, supplemented with 5% FBS, penicillin (200 U ml<sup>-1</sup>), streptomycin (200  $\mu$ g ml<sup>-1</sup>) and 4 mM L-glutamine and passaged using Trypsin-EDTA (0.05%). Both SK-BR-3 and LNCaP cells were grown in RPMI 1640 medium, supplemented with 20% FBS, penicillin (200 U ml<sup>-1</sup>), streptomycin (200  $\mu$ g ml<sup>-1</sup>) and 2 mM L-glutamine and passaged using TrypLE Select Enzyme.

The murine A20 cell line was stored in liquid nitrogen and grown in RPMI 1640, supplemented with 10% FBS, penicillin (200 U ml<sup>-1</sup>), streptomycin (200  $\mu$ g ml<sup>-1</sup>), 200  $\mu$ M L-glutamine, 25 mM D-glucose, 10 mM HEPES and 1 mM sodium pyruvate.

The cells were cultivated in T-75 cell culture flasks at 37 °C and 5% CO<sub>2</sub>. The cells were maintained according to the instructions from DSMZ and used for flow cytometric and cell-based assays up to a maximum of 15 passages. This cell line was used as a murine B-cell lymphoma model.

For in vivo models, NALM-6 tumour cells were lentivirally transduced with a pCDH-EF1a-eFly-eGFP plasmid as previously described<sup>49,50</sup>. Short tandem repeat profiling was used to verify the origin of this cell line.

For flow cytometric experiments, the following antibodies were used as purchased.

Antibody	Company
CD8a monoclonal antibody (RPA-T8)	Thermo Fisher
CD69 monoclonal antibody (FN50)	Thermo Fisher
PerCP-eFluor 710 CD4 monoclonal antibody (RPA-T4), PE	Thermo Fisher
CD25 monoclonal antibody (PC61.5), PE-eFluor 610	Thermo Fisher

#### Folding of DNA origami objects (chassis)

The reaction mixtures contained scaffold DNA at a concentration of 50 nM and oligonucleotide strands at 200 nM each. The reaction buffer included 5 mM Tris, 1 mM EDTA, 5 mM NaCl (pH 8) and 20 mM MgCl<sub>2</sub>. The reaction mixtures were subjected to a thermal annealing ramp using Tetrad (MJ Research, now Bio-Rad) thermal cycling devices. Oligonucleotides were obtained from IDT. DNA scaffolds were produced in-house according to another work<sup>51</sup>. The table below shows the folding ramps used to assemble the objects described in the manuscript.

Object	Scaffold type	Highest temperature for 30 min (°C)	Ramp	Incubation temperature (°C)
Medium brick	8,064	65	55 to 50; 1°C/1h	25
Small brick	1,033	65	56 to 54; 1°C/2h	25

#### Gel electrophoresis of PTEs

Folded DNA nanostructures were electrophoresed on 1.5% to 3.5% agarose gels containing  $0.5 \times$  TBE and MgCl<sub>2</sub> at different concentrations for around 2 h at 70 V bias voltage in a gel box immersed in a water bath, unless specified otherwise. The electrophoresed agarose gels were scanned using a Typhoon FLA 9500 laser scanner (GE Healthcare) at 50 µm per pixel resolution. The resulting 16-bit TIFF images were analysed using ImageJ v. 1.440. For each lane that contained the sample, a cross-sectional intensity profile was calculated by averaging over grey-scale values within a 50-pixel-wide box. The peak intensities of the monomers and higher-ordered bands were determined in the target band. These intensity values have been used for further analysis.

#### Purification, enrichment and in vitro stabilization of PTEs

After the folding reaction, all the reaction products were purified using PEG precipitation<sup>52</sup>. For concentrating DNA origami objects (chassis or PTEs), PEG precipitation or ultrafiltration was used. All the procedures were performed as previously described<sup>53</sup>. Concentrations of DNA origami objects were analysed with a Nanodrop 8000 instrument (Thermo Fisher). Before using the objects in cell culture assays, the objects were stabilized for use in low-ionic-strength buffers and the presence of nucleases. To that end, we used the protocol from another work<sup>22</sup> and coated all our structures with K10-PEG oligolysine, purchased from Alamanda Polymers.

# Antibody–DNA conjugation, antibody digestion and attachment to DNA origami objects

Conjugation of full-sized IgG antibodies: oligonucleotides modified with 5'- or a 3'-thiol modification were purchased, HPLC purified and dried (Biomers). The oligos were dissolved in PBS (100 mM NaPi, 150 mM NaCl, pH 7.2) with 5 mM TCEP and incubated for 1 h at RT. After purification, 10 nmol of the reduced thiol oligo was mixed with 10 equivalents of sulfosuccinimidyl 4-(*N*-maleimidomethyl) cyclohexane-1-carboxylate (Sulfo-SMCC; dissolved in ddH<sub>2</sub>O) (Thermo Fisher) for 15 min. After purification, including buffer change to PBS (pH 8), 100 µg of antibody in PBS (pH 8) was added. The conjugate was subsequently purified by ion exchange chromatography (proFIRE, Dynamic Biosensors) using a NaCl gradient of 150–1,000 mM in PBS (pH 7.2). Purified oligo–antibody conjugates were analysed by sodium dodecyl sulfate–polyacrylamide gel electrophoresis and agarose gel electrophoresis.

Preparation of Fab fragments: full-sized IgGs were digested and purified using the Pierce Fab Preparation Kit (Thermo Fisher, 44985) according to the supplier's manual. Briefly, IgGs were incubated with papain beads and purified from Fc fragments using Protein-A affinity beads. Fab generation was checked using sodium dodecyl sulfate– polyacrylamide gel electrophoresis.

Conjugation of Fab fragments: to avoid an orientation where the F(ab)'s paratope points towards the DNA origami chassis, we relied on a site-specific conjugation method similar to another work<sup>54</sup>. In brief, maleimide-modified DNA was purchased from Biomers or prepared by mixing amine-modified DNA with an SMCC crosslinker and subsequent ultracentrifugation (10k filters, Amicon). Fab fragments, which were produced using the papaine digestion of IgG (Pierce Fab Preparation Kit, Thermo Fisher, 44985) were reduced with 5 mM TCEP for 30 min. Excess TCEP was removed using ultracentrifugation (10k filters, Amicon), and mixed with maleimide-modified DNA strands in 50 mM

HEPES with 200 mM NaCl at pH 6.7. Reactions were performed overnight at room temperature. This conjugation technique results in single- and double-labelled Fab fragments. Single-labelled F(ab)–DNA conjugates were purified using ion exchange chromatography.

Conjugation of scFv: maleimide-modified DNA was purchased from Biomers or prepared by mixing amine-modified DNA with SMCC crosslinker. scFv contained a free N-terminal cysteine and were expressed and purchased from Genscript or Icosagen. The reaction was carried out as described for the Fab fragments.

Attachment of antibody–DNA conjugates to DNA origami objects: antibody–DNA conjugates and DNA origami objects with the corresponding binding sites were incubated in equimolar ratios for 1 h at 37 °C.

#### **Negative-stain TEM**

**Preparation, acquisition and data processing.** Purified reaction products were adsorbed on glow-discharged Cu400 TEM grids (Science Services) and stained using a 2% aqueous uranyl formate solution containing 25 mM sodium hydroxide. The samples were incubated for 30 s at 20-25 mM Mg<sup>2+</sup>. Magnifications between ×10,000 to ×30,000 were used for acquiring the data.

Imaging was performed on different microscopes.

Microscope	Operating voltage (kV)	Camera	Objects
Philips CM 100	100	AMT 4 megapixel charge-coupled device camera	Bricks S, M
FEI Tecnai 120	120	TemCam F416 (4k×4k)	Bricks S, M

TEM micrographs used in the figures were high-pass filtered to remove long-range staining gradients and the contrast was auto-levelled (Adobe Photoshop CS6).

#### Fluorescence microscopy

Fluorescence microscopy experiments were obtained on a Thermo Fisher CX7 confocal microscope with an on-stage incubator. Incubation conditions for time-resolved measurements were identical to the cell culture conditions used for the respective cell lines. The samples were incubated on 96-well ibidi plates (89626).

#### Cell-based assays

**Cell surface binding experiments.** For flow cytometric experiments, the cells were grown to a cell density of  $1.5-2.0 \times 10^6$  cells ml<sup>-1</sup> in T-75 cell culture flasks. The cells were centrifuged for 5 min at  $160 \times g$  and washed with ice-cold PBS, twice. All the flow cytometric experiments were executed at a cell density of  $2 \times 10^7$  cells ml<sup>-1</sup> in PBS or medium. The sample (chassis or PTE) was added at a final concentration of 1 nM and incubated for the different time points. Before flow cytometric analysis, the cells were centrifuged for 5 min at  $500 \times g$  and resuspended to a final coll concentration of  $2 \times 10^6$  cells ml<sup>-1</sup> in PBS. Flow cytometric analysis was performed on a Cytoflex (Beckman Coulter) or Attune Nxt (Thermo Fisher), measuring the fluorescent intensity by excitation at 640 nm and band-pass filter detection at 660/20 nm. The single cells were gated based on the forward scatter versus side scatter. For each measurement, fluorescent intensities of 50,000 individual cells were analysed with an in-house MATLAB (R2021b) script.

**T-cell activation assay.** The expression of interleukin 2, as an indicator for T-cell activation, was analysed using T-Cell Activation Bioassays (Promega, J1655)<sup>35</sup>. The experiment was performed according to the supplier's instructions. Briefly, CD19-expressing target cells (NALM-6) were added to 96-well microtiter plates at a final concentration of  $5 \times 10^5$  cells ml<sup>-1</sup>. Then, a serial dilution of different samples (in RPMI 1640 medium) was added. In the end, the genetically modified TCR/CD3

effector cells were added at a final concentration of  $1.3 \times 10^6$  cells ml<sup>-1</sup>. The reaction mixture was incubated for 6 h at 37 °C and 5% CO<sub>2</sub>. The genetically modified effector cells (Jurkat-NFAT) intercellularly express a luciferase if the interleukin 2 promoter is activated. By the addition of the Bio-Glo reagent, which includes a substrate for the luciferase, the luminescence signal is a direct proportional signal for the activation of the TCR/CD3 effector cells, which was analysed in a microtiter plate reader (Clariostar Plus, BMG). Data were normalized and the background signal was corrected.

**Internalization assay.** The internalization assay was performed according to another work<sup>55</sup>. Each DNA origami brick (chassis) carried a fluorescence internalization probe (FIP) comprising a protruding sequence with a terminal-attached Cy5 dye (Supplementary Fig. 4, right inset). This FIP can be quenched using a quenching strand with a complementary sequence to the FIP and an attached Black Hole Quencher-2. On hybridization, Black Hole Quencher-2 quenches the fluorescence of the FIP. Since the quencher strand can only reach the chassis on the cell surface, and not the internalized chassis, the amount of quenched signal is proportional to the amount of surface-exposed chassis.

NALM-6 cells ( $1 \times 10^7$  cells ml<sup>-1</sup>) were incubated with 1 nM chassis for 1 h at 37 °C in the cell culture medium and then washed to remove the excess chassis. After washing, the cells were resuspended in the cell culture medium and incubated at 37 °C. At each time point, a measurement consists of taking a sample and incubating it for 10 min at 4 °C on ice. Half of the sample is incubated without a quencher strand and the other half is incubated with a 100 nM quencher strand. Both samples were incubated for 10 min on ice to allow for quencher hybridization (if added) and to stop internalization. The fraction of the internalized chassis with two anti-CD19 antibodies was calculated from the median fluorescence F as follows.

Fraction internalized =

 $\frac{F_{\text{with antibodies, with quencher}} - F_{\text{without antibodies, with quencher}}}{F_{\text{with antibodies, without quencher}} - F_{\text{without antibodies, without quencher}}}$ 

**Cytotoxic T-cell-killing assay of liquid tumour cells.** Preparation of target cells. Target cells (for example, NALM-6) were fluorescently stained using CellTrace CSFE Cell Proliferation Kit (Thermo Fisher, C34554) according to the manufacturer's protocol.

**Preparation of PBMC (effector).** For the cytotoxic T-cell-killing assay, we used frozen human PBMC (STEMCELL, 70025.1 or CTL, CTL-UP1). We handled the PBMC cells according to the supplier's instruction.

**Assay.** Here  $2 \times 10^{5}$  CSFE-stained target cells per millilitre were incubated with  $1 \times 10^{6}$  PBMC ml<sup>-1</sup> in the cell culture medium at 37 °C (5% CO<sub>2</sub>) and the PTEs in different concentrations or without an additional recruiter. Cell fluorescence and scattering intensity were determined using an Attune Nxt flow cytometry with a Cytkick Max autosampler (Thermo Fisher). Residual contaminations, such as salts or endotoxins may cause the non-specific lysis of target cells. These effects are particularly pronounced at high PTE concentrations and depend on the cell type.

**Cytotoxic T-cell-killing assay of solid tumour cells.** Preparation of target cells: target cells (for example, MCF-7) were seeded 24 h in advance to from a confluent monolayer.

Preparation of PBMC (effector): for the cytotoxic T-cell-killing assay, we used frozen human PBMC (STEMCELL, 70025.1 or CTL, CTL-UP1). We handled the PBMC cells according to the supplier's instruction.

Assay: confluent target cells were incubated with PBMC (E:T 2:1 or 4:1) in cell culture medium at 37 °C (5%  $CO_2$ ) and the PTEs in different concentrations or without an additional recruiter. The fraction of alive target cells was determined by quantifying the amount of ATP after 48 h

via Bio-Glo Cell Titer System (Promega) of samples with and without PTEs. Residual contaminations, such as salts or endotoxins, may cause the non-specific lysis of target cells. These effects are particularly pronounced at high PTE concentrations and depend on the cell type.

#### In vivo animal models

**Endotoxin determination of in vivo ready constructs.** The endotoxin concentration was measured with a Charles River nexgen-PTS150 V10.2.3 instrument. We used cartridges with a range between 5.00 and 0.05 EU ml<sup>-1</sup>. The samples were diluted 25-fold to fit into the sensitive range of the cartridges. The endotoxin threshold level for mouse studies was set to 36 EU ml<sup>-1</sup>. This value is in accordance with the specifications given by the FDA<sup>42</sup>.

Sample	EU measured	EU in sample (25-fold)	QC passed
1-0 (origami)	0.068	1.7	Yes
0-2 (origami)	0.068	1.7	Yes
1-2 (origami)	0.083	2.0	Yes

**In vivo mouse experiments.** Approval for all the animal experiments was granted by the local regulatory authorities (Regierung von Oberbayern).

NOD.Cg-Prkdc<sup>scid</sup> Il2rg<sup>tm1WjI/Szl</sup> (NSG) were purchased from Janvier or bred in-house. NSG mice carry a mutation in the Pkrdc DNA repair gene, associated with severe combined immunodeficiency, leading to a T- and B-cell deficiency. The complete null allele of the IL-2 receptor gamma chain (Il2rg) abrogates the cytokine signalling of critical homoeostatic cytokines such as IL-2, IL-4, IL-7, IL-9 and IL-15, preventing the development of functional NK cells. Finally, the NOD background further compromises the innate branch of the immune system (reduced functionality of dendritic cells and macrophages). In general, these highly immunodeficient mice support stable and reproducible engraftment of human tumour and T cells in mice and are currently regarded as the state-of-the-art model for human xenograft models in mice<sup>56–58</sup>.

All the animals were housed in specific pathogen-free facilities.

BLI and fluorescence imaging were carried out using the in vivo imaging platform Lumina X5 (IVIS, PerkinElmer), as previously described<sup>59</sup>. In brief, mice were anaesthetized with a 1.5–2.5% isoflurane-oxygen mixture for all the live-imaging procedures. For BLI, the substrate (Xenolight D-luciferin potassium salt, PerkinElmer) was intraperitoneally injected according to the manufacturer's instructions. For organ analysis, the background fluorescence of each organ was subtracted.

NALM-6-luc $^+$ -GFP $^+$  xenograft models were established by i.v. tail-vein injection.

No data points (mice) were excluded in the animal studies.

#### Flow cytometry for in vivo experiments

Flow cytometric data were generated using a BD LSRFortessa II, a Beckman Coulter CytoFLEX LX or a Thermo Fisher Attune Nxt with an autoloader. Flow cytometric analysis of organs was conducted as previously reported<sup>60</sup>. Single-cell suspensions of harvested organs were stained with human anti-CD3 BV711 (clone, OKT3), anti-CD4 PerCP-Cy5.5 (clone, OKT4), anti-CD8 PE (clone, HIT8a), anti-CD19 BV786 (clone, HIB19) and anti-CD69 PE-Cy7 (clone, FN50) or mouse anti-CD45 pacific blue (clone, 30-F11) antibodies (Biolegend). Fixable viability dye (eFluor 780, eBioscience) was used to exclude the dead cells. The maximum tumour burden permitted by the local regulatory authorities was not exceeded.

**Software and statistical analyses.** Flow cytometric data were analysed using FlowJo v. 10.3 to v. 10.8.1 software. Quantifying bioluminescence and fluorescence intensities was done using Living Image 4.4

(PerkinElmer). Statistical analysis was carried out with GraphPad Prism v. 9.4.0. Power calculations (for in vivo experiments) were carried using G\*Power 3.1 with given alpha, power and effect size.

#### Statistics and reproducibility

Statistics and reproducibility are stated in the legend of the figures, for example, biological or technical replicates. Agarose gel images shown in the figures are representative examples of experiments that yielded the same or similar results. For TEM analysis: the total number of similarly conducted TEM analysis of samples prepared following the same protocols varied between experiments—in case of replicates, reproducibility was observed.

#### **Reporting summary**

Further information on research design is available in the Nature Portfolio Reporting Summary linked to this article.

#### **Data availability**

All data that support the findings of this study are available within the Article and its Supplementary Information, and available from the corresponding authors on request. Source data are provided with this paper.

#### References

- Sandhöfer, N. et al. Dual PI3K/mTOR inhibition shows antileukemic activity in MLL-rearranged acute myeloid leukemia. *Leukemia* 29, 828–838 (2015).
- Karches, C. H. et al. Bispecific antibodies enable synthetic agonistic receptor-transduced T cells for tumor immunotherapy bispecific antibodies enable transduced T cells. *Clin. Cancer Res.* 25, 5890–5900 (2019).
- 51. Engelhardt, F. A. et al. Custom-size, functional, and durable DNA origami with design-specific scaffolds. *ACS Nano* **13**, 5015–5027 (2019).
- 52. Stahl, E., Martin, T. G., Praetorius, F. & Dietz, H. Facile and scalable preparation of pure and dense DNA origami solutions. *Angew. Chem.* **53**, 12735–12740 (2014).
- 53. Wagenbauer, K. F. et al. How we make DNA origami. *ChemBioChem* **18**, 1873–1885 (2017).
- 54. Humphreys, D. P. et al. Alternative antibody Fab' fragment PEGylation strategies: combination of strong reducing agents, disruption of the interchain disulphide bond and disulphide engineering. *Protein Eng. Des. Sel.* **20**, 227–234 (2007).
- 55. Liu, H. & Johnston, A. P. A programmable sensor to probe the internalization of proteins and nanoparticles in live cells. *Angew. Chem.* **125**, 5856–5860 (2013).
- McDermott, S. P., Eppert, K., Lechman, E. R., Doedens, M. & Dick, J. E. Comparison of human cord blood engraftment between immunocompromised mouse strains. *Blood* **116**, 193–200 (2010).
- Shultz, L. D. et al. Human lymphoid and myeloid cell development in NOD/LtSz-scid IL2R gamma null mice engrafted with mobilized human hemopoietic stem cells. *J. Immunol.* **174**, 6477–6489 (2005).
- Lepus, C. M. et al. Comparison of human fetal liver, umbilical cord blood, and adult blood hematopoietic stem cell engraftment in NOD-scid/γc<sup>-/-</sup>, Balb/c-Rag1<sup>-/-</sup>γc<sup>-/-</sup>, and C.B-17-scid/bg immunodeficient mice. Hum. Immunol. **70**, 790–802 (2009).
- 59. Benmebarek, M.-R. et al. A modular and controllable T cell therapy platform for acute myeloid leukemia. *Leukemia* **35**, 2243–2257 (2021).
- 60. Lesch, S. et al. T cells armed with CXC chemokine receptor type 6 enhance adoptive cell therapy for pancreatic tumours. *Nat. Biomed. Eng.* **5**, 1246–1260 (2021).

#### Acknowledgements

This work was supported by a European Research Council Consolidator Grant to H.D. (grant agreement 724261), the Deutsche Forschungsgemeinschaft (DFG) through grants provided within the Gottfried Wilhelm Leibniz Program (to H.D.), the Medical Valley Award, the M4 Award, a GO-Bio initial award (Federal Ministry of Education and Research (BMBF) of Germany) and funding by the ForTra gGmbH für Forschungstransfer der Else Kröner-Fresenius Stiftung to K.F.W., J.J.F., B.K. and H.D. The work has received support from the Max Planck School Matter to Life (to H.D.) jointly financed by the Federal Ministry of Education and Research (BMBF) of Germany and the Max Planck Society. K.F.W. is grateful for support from the Bavarian Academy of Science. This study was further supported by the international doctoral program the Förderprogramm für Forschung und Lehre der Medizinischen Fakultät der LMU (to A.G.; grant no. 1139), the DFG (to A.G. via grant no. GO 3823/1-1 and to S.K. via grant no. KO5055/3-1), 'i-Target: immunotargeting of cancer' (funded by the Elite Network of Bavaria; to S.K.), Melanoma Research Alliance (grant no. 409510 to S.K.), Marie Skłodowska-Curie Training Network for Optimizing Adoptive T Cell Therapy of Cancer (funded by the Horizon 2020 programme of the European Union; grant 955575 to S.K.), Else Kröner-Fresenius-Stiftung (to A.G. and S.K.), German Cancer Aid (AvantCAR.de to S.K.), Ernst Jung Stiftung (to S.K.), the Wilhelm-Sander Stiftung (to S.K.), Institutional Strategy LMUexcellent of LMU Munich (within the framework of the German Excellence Initiative; to S.K.), Bundesministerium für Bildung und Forschung (to S.K.), European Research Council (Starting Grant 756017 and Proof of Concept Grant 101100460 to S.K.), DFG (KO5055-2-1 and 510821390 to S.K.), by the SFB-TRR 338/1 2021-452881907 (to S.K.), Fritz-Bender Foundation (to S.K.), Deutsche José Carreras Leukämie Stiftung (to S.K.) and Hector Foundation (to S.K.). Cytometry data were obtained in the Core Facility Flow Cytometry of the University Hospital, LMU Munich, using a BD LSRFortessa II or Beckman Coulter CytoFLEX LX. The in vivo imaging device was funded by the DFG (German Research Foundation)—INST 409/231-1.

#### **Author contributions**

K.F.W., B.K., J.J.F. and H.D. developed the concept. K.F.W., A.G., B.K., S.K., J.J.F. and H.D. designed the experiments. K.F.W., N.P., B.K.,

V.K., C.F., P.S. and J.J.F. performed the experiments in vitro. A.G., R.G., E.C. and C.A.T. performed the animal experiments. K.F.W., A.G., B.K., J.J.F. and H.D. wrote the manuscript. K.F.W., B.K., S.K., J.J.F. and H.D. supervised the project and acquired the funding. All authors critically read and approved the manuscript.

#### **Competing interests**

The Technical University of Munich has filed several patents listing K.F.W., J.J.F., B.K. and H.D. as the inventors. K.F.W., J.J.F., B.K. and H.D. are co-founders of Plectonic Biotech GmbH. S.K. has received honoraria from TCR2 Inc., Novartis, Miltenvi Biomedicines, BMS and GSK, S.K. and A.G. are inventors of several patents in the field of immuno-oncology (approved, S.K., PCT/EP2013/051351, PCT/EP2016/064195, PCT/EP2016/074644; pending, S.K., PCT/EP2020/056086; submitted, S.K. and A.G., EP21191376, EP 23 154 047.78; countries, Europe, USA, Canada). S.K. received license fees from TCR2 Inc. and Carina Biotech. S.K. received research support from Plectonic GmbH for parts of the work contained herein. S.K. received research support from TCR2 Inc. and Arcus Bioscience for work unrelated to the manuscript. A.G. received research support from Tabby Therapeutics for work unrelated to the manuscript. The remaining authors declare no competing interests.

#### **Additional information**

**Supplementary information** The online version contains supplementary material available at https://doi.org/10.1038/s41565-023-01471-7.

**Correspondence and requests for materials** should be addressed to Sebastian Kobold, Jonas J. Funke or Hendrik Dietz.

**Peer review information** *Nature Nanotechnology* thanks Eva Sevcsik and the other, anonymous, reviewer(s) for their contribution to the peer review of this work.

**Reprints and permissions information** is available at www.nature.com/reprints.

# nature portfolio

Corresponding author(s): instead of author names.

Last updated by author(s): Jun 9, 2023

# **Reporting Summary**

Nature Portfolio wishes to improve the reproducibility of the work that we publish. This form provides structure for consistency and transparency in reporting. For further information on Nature Portfolio policies, see our <u>Editorial Policies</u> and the <u>Editorial Policy Checklist</u>.

### **Statistics**

For	all st	atistical analyses, confirm that the following items are present in the figure legend, table legend, main text, or Methods section.
n/a	Cor	firmed
	$\square$	The exact sample size (n) for each experimental group/condition, given as a discrete number and unit of measurement
	$\square$	A statement on whether measurements were taken from distinct samples or whether the same sample was measured repeatedly
		The statistical test(s) used AND whether they are one- or two-sided Only common tests should be described solely by name; describe more complex techniques in the Methods section.
$\boxtimes$		A description of all covariates tested
$\boxtimes$		A description of any assumptions or corrections, such as tests of normality and adjustment for multiple comparisons
	$\boxtimes$	A full description of the statistical parameters including central tendency (e.g. means) or other basic estimates (e.g. regression coefficient) AND variation (e.g. standard deviation) or associated estimates of uncertainty (e.g. confidence intervals)
		For null hypothesis testing, the test statistic (e.g. F, t, r) with confidence intervals, effect sizes, degrees of freedom and P value noted Give P values as exact values whenever suitable.
$\ge$		For Bayesian analysis, information on the choice of priors and Markov chain Monte Carlo settings
$\boxtimes$		For hierarchical and complex designs, identification of the appropriate level for tests and full reporting of outcomes
$\boxtimes$		Estimates of effect sizes (e.g. Cohen's d, Pearson's r), indicating how they were calculated
		Our web collection on <u>statistics for biologists</u> contains articles on many of the points above.

## Software and code

 Policy information about availability of computer code

 Data collection
 Data were collected on different instruments. List of instruments is included in the methods section in the S.I.

 Nanodrop 8000 Thermo Fisher, Attune Nxt Thermo Fisher, Typhoon GE, Profire Dynamic Biosensors, Platereader Clariostar BMG Labtech, Charles River NexGen PTS, CX7 Thermo Fisher, BD LSR Fortessa II, a Beckman Coulter CytoFLEX LX, Living Image 4.4 (PerkinElmer, USA)

 Data analysis
 Data analysis is described in detail in the S.I. Data were analyzed with the following software Matlab 2021b, Igor Pro 7.06, Office 365, Fiji 2.1.0, FlowJo V10.3 to V10.8.1 software, GraphPad Prism V.9.4.0 (San Diego, CA, USA).

For manuscripts utilizing custom algorithms or software that are central to the research but not yet described in published literature, software must be made available to editors and reviewers. We strongly encourage code deposition in a community repository (e.g. GitHub). See the Nature Portfolio guidelines for submitting code & software for further information.

Policy information about <u>availability of data</u>

All manuscripts must include a <u>data availability statement</u>. This statement should provide the following information, where applicable:

- Accession codes, unique identifiers, or web links for publicly available datasets
- A description of any restrictions on data availability
- For clinical datasets or third party data, please ensure that the statement adheres to our policy

All data and source data that support the findings of this study are available within the paper and its Supplementary Information, and from the corresponding author on request.

## Research involving human participants, their data, or biological material

Policy information about studies with <u>human participants or human data</u>. See also policy information about <u>sex, gender (identity/presentation)</u>, <u>and sexual orientation</u> and <u>race, ethnicity and racism</u>.

Reporting on sex and gender	N.A.
Reporting on race, ethnicity, or other socially relevant groupings	N.A.
Population characteristics	N.A.
Recruitment	N.A.
Ethics oversight	N.A.

Note that full information on the approval of the study protocol must also be provided in the manuscript.

# Field-specific reporting

Please select the one below that is the best fit for your research. If you are not sure, read the appropriate sections before making your selection.

Life sciences Behavioural & social sciences Ecological, evolutionary & environmental sciences

For a reference copy of the document with all sections, see <u>nature.com/documents/nr-reporting-summary-flat.pdf</u>

# Life sciences study design

All studies must disclose on these points even when the disclosure is negative.

Sample size	For animal experiments sample size was determined using the software G*Power 3.1 with given alpha, power and effect size	
Data exclusions	No data exclusion was done.	
Replication	For TEM and gel electrophoretic analysis replicates were performed in different arangements. Images in the manuscript are examples of such experiments. For in-vitro T-cell mediated lysis experiments in general biological triplicates were done if not otherwise stated and explainded. A20 cell line experiments were done once with Splenocytes due to ethical reasons. T-cell activation experiments were done in technical dublicates. For in vivo experiments, group size and number of independent repetitions are indicated. In vivo experiments were performed once (Fig. 4,	
	Supplementary Figure 19, 20) or three-times (Fig. 3).	
Randomization	Randomization was performed for all animal experiments. If tumor sizes showed strong variations, mice were stratified to ensure similiar tumor burden across all treatment groups.	
Blinding	All animal experiments were performed blinded.	

# Reporting for specific materials, systems and methods

We require information from authors about some types of materials, experimental systems and methods used in many studies. Here, indicate whether each material, system or method listed is relevant to your study. If you are not sure if a list item applies to your research, read the appropriate section before selecting a response.

#### Materials & experimental systems

n/a	Involved in the study
	Antibodies
	Eukaryotic cell lines
$\boxtimes$	Palaeontology and archaeology
	Animals and other organisms
$\boxtimes$	Clinical data
$\boxtimes$	Dual use research of concern
$\boxtimes$	Plants

#### Methods

n/a Involved in the study

Flow cytometry

MRI-based neuroimaging

## Antibodies

Antibodies used	Antibody Clon Vendor ID CD3 UTCH1 ThermoFisher 16-0038-85 CD3 OKT3 Biozol BE0001-2 CD3 OKT3 losagen Custom - scFv CD19 HIB19 ThermoFisher 16-0199-85 CD123 6H6 ThermoFisher 14-1239-37 CD22 eBio4KB128 (4KB128) ThermoFisher 14-0229-82 CD33 WM-53 ThermoFisher 14-0338-37 CD28 CD28.2 ThermoFisher 16-0289-85 CD137 5G11 ThermoFisher 16-0289-85 CD137 5G11 ThermoFisher 14-9056-82 muCD3 F(ab)2 145-2C11 Biozol 711721J1 EpcAM VU-109 ThermoFisher MA5-12153 PSMA GCP-05 ThermoFisher MA1-10335 CD8a Monoclonal Antibody (RPA-T8) ThermoFisher CD69 Monoclonal Antibody (RPA-T8) ThermoFisher CD59 Monoclonal Antibody (RPA-T8), PE-eFluor™ 710 ThermoFisher CD55 Monoclonal Antibody (PC61.5), PE-eFluor™ 610 ThermoFisher Flow cytometry: TrueStain FcX™ (BioLegend, USA, catalog #422302) fixable viability dye eFluor™ 780 (eBioscience, USA, catalog #65-0865-18) anti-murine CD45 - PacBlue (Biolegend, clone: 30-F11, catalog #103126) anti-human CD69 PE-Cy7 ( human – clone: FN50, catalog #310912, Biolegend, USA) anti-CD3 (Biolegend, BV711 - clone: UCHT1, catalog 300464) anti-CD4 (Biolegend, human, PECP-Cy5.5 - clone: OXT4, catalog #317428) anti-CD8 (Biolegend, human, PE – clone HIT8a, catalog #309098)
Validation	CD3 UTCH1 ThermoFisher 16-0038-85 https://www.thermofisher.com/antibody/product/CD3-Antibody-clone-UCHT1- Monoclonal/16-0038-85 CD3 OKT3 Biozol BE0001-2 https://www.biozol.de/InVivoMab-anti-human-CD3-Clone-OKT-3-Mouse-Monoclonal/BXC-BE0001-2-1MG CD3 OKT3 Biozol BE0001-2 https://www.thermofisher.com/antibody/product/CD19-Antibody-clone-HIB19- Monoclonal/16-0199-85 CD123 6H6 ThermoFisher 16-0199-85 https://www.thermofisher.com/antibody/product/CD123-Antibody-clone-HIB19- Monoclonal/14-1239-37 CD22 eBio4KB128 (4KB128) ThermoFisher 14-0229-82 https://www.thermofisher.com/antibody/product/CD123-Antibody-clone-6H6- Monoclonal/14-1239-37 CD22 eBio4KB128 (4KB128) ThermoFisher 14-0229-82 https://www.thermofisher.com/antibody/product/CD23-Antibody-clone- eBio4KB128-4KB128-Monoclonal/14-0229-82 CD33 WM-53 ThermoFisher 14-0338-37 https://www.thermofisher.com/antibody/product/CD23-Antibody-clone- CD28 CD28 2. ThermoFisher 14-038-37 https://www.thermofisher.com/antibody/product/CD23-Antibody-clone-CD28-2- Monoclonal/14-0338-37 CD127 5G11 ThermoFisher 14-9056-82 https://www.thermofisher.com/antibody/product/CD137-Ligand-4-1B8-Ligand-Antibody- clone-5G11-Monoclonal/14-9056-82 mUCD19 1D3 Biozol BE0150 https://www.biozol.de/InVivoMab-anti-mouse-CD19-Clone-1D3-Rat-Monoclonal/BXC-BE0150-5MG muCD3 F(ab)2 145-2C11 Biozol 71172111 https://www.thermofisher.com/antibody/product/CD137-Ligand-4-1B8-Ligand-Antibody- Clone-5G11-Monoclonal/X4-9056-82 mUCD19 1D3 Biozol BE0150 https://www.thermofisher.com/antibody/product/EpCAM-Antibody-clone-VU-1D9- Monoclonal/MA5-12153 PSMA GCP-05 ThermoFisher MA5-12153 https://www.thermofisher.com/antibody/product/PSMA-Antibody-clone-GCP-05- Monoclonal/MA1-10335 CD8a Monoclonal Antibody (RPA-T8) https://www.thermofisher.com/antibody/product/CD8a-Antibody-clone-RPA-T8- Monoclonal/MA1-10335 CD69 Monoclonal Antibody (RPA-T8) https://www.thermofisher.com/antibody/product/CD4-Antibody-clone-RPA-T8- Monoclonal/A17-0088-42 CD69 Monoclonal Antibody (RPA-T4), PE https://www.thermofisher.com/antibody/product/CD4-An

MA1-81104

CD25 Monoclonal Antibody (PC61.5), PE-eFluor<sup>™</sup> 610 https://www.thermofisher.com/antibody/product/CD25-Antibody-clone-PC61-5-Monoclonal/61-0251-80

TrueStain FcX™ (https://www.biolegend.com/en-us/products/human-trustain-fcx-fc-receptor-blocking-solution-6462? GroupID=BLG2181) fixable viability dye eFluor™ 780 (eBioscience, USA, catalog #65-0865-18)

anti-murine CD45 - PacBlue (https://www.biolegend.com/en-us/products/pacific-blue-anti-mouse-cd45-antibody-19250? GroupID=GROUP20)

anti-human CD69 PE-Cy7 (https://www.biolegend.com/en-us/products/pe-cyanine7-anti-human-cd69-antibody-1918)

anti-CD3 (ttps://www.biolegend.com/en-us/products/brilliant-violet-711-anti-human-cd3-antibody-12047) anti-human CD4 PerCP-Cy5.5 (https://www.biolegend.com/ja-jp/products/percp-cyanine5-5-anti-human-cd4-antibody-5011) anti-human CD8 PE (https://www.biolegend.com/en-ie/products/pe-anti-human-cd8a-antibody-762)

## Eukaryotic cell lines

Policy information about <u>cell lines</u>	and Sex and Gender in Research
Cell line source(s)	Jurkat https://www.dsmz.de/collection/catalogue/details/culture/ACC-282 NALM-6 https://www.dsmz.de/collection/catalogue/details/culture/ACC-128 MCF-7 https://www.dsmz.de/collection/catalogue/details/culture/ACC-115 Molm-13 https://www.dsmz.de/collection/catalogue/details/culture/ACC-554 LNCaP https://cls.shop/LNCaP/300265 NALM-6 (ATCC, USA)
Authentication	Cell line Authentification was done by vendor and via PCR by Eurofins Genomics. STR DNA profiling of human cell lines (NALM-6).
Mycoplasma contamination	All cell lines were tested negative for mycoplasma contamination.
Commonly misidentified lines (See <u>ICLAC</u> register)	No commonly misidentified lines were used.

## Animals and other research organisms

Policy information about studies involving animals; <u>ARRIVE guidelines</u> recommended for reporting animal research, and <u>Sex and Gender in</u> <u>Research</u>

Laboratory animals	NSG (NOD.Cg-PrkdcSCIDII2rgtm1Wjl/SzJ mice, Charles River, Janiver or breeding in the local animal facility "Zentrale Versuchstierhaltung Innenstadt", male or female, 4-12 weeks and >16g body weight at the start of the experiments). Mice were held in facilities with a 12-hours dark/light cycle including a 30 min twilight phase of 30 minutes at noise levels below 50 dBA. Air velocity was held below 0.2 m/s. Air humidity in the facilities were between 45 – 60% and average temperature was held between 20 to 22 °C.
Wild animals	No wild animals were used in the study.
Reporting on sex	Findings are applicable to both male and female mice.
Field-collected samples	No field collected samples were used in the study.
Ethics oversight	Local regulatory agency (Regierung von Oberbayern)

Note that full information on the approval of the study protocol must also be provided in the manuscript.

## Flow Cytometry

#### Plots

Confirm that:

The axis labels state the marker and fluorochrome used (e.g. CD4-FITC).

The axis scales are clearly visible. Include numbers along axes only for bottom left plot of group (a 'group' is an analysis of identical markers).

 $\square$  All plots are contour plots with outliers or pseudocolor plots.

A numerical value for number of cells or percentage (with statistics) is provided.

#### Methodology

Sample preparation

Cell surface binding experiments

Sample preparation	cells were centrifuged for 5 min at 160 rcf and washed with ice-cold PBS, twice. All flow cytometry experiments were executed at a cell density of 2*10^7 cells ml-1 in PBS or medium. The sample (Chassis or PTE) was added at a final concentration of 1 nM and incubated for the different time points. Before flow cytometry analysis, the cells were centrifuged for 5 min at 500 rcf and resuspended to a final cell concentration of 2*10^6 cell ml-1 in PBS. Flow cytometric analysis was performed on a Cytoflex (Beckman Coulter) or Attune Nxt (ThermoFisher), measuring the fluorescent intensity by excitation at 640 nm and a bandpass detection filter at 660/20nm. The single cells were gated based on the forward scatter versus the side scatter. Cytotoxic T-cell killing assay of liquid tumor cells: Assay: 2x10^5 CSFE stained target cells per ml were incubated with 1x10^6 PBMC/ml in the cell-culture medium at 37°C (5% CO2) and the programmable T-cell engagers in different concentrations or without an additional recruiter. Cell fluorescence and scattering intensity was determined using an Attune Nxt flow cytometry with a Cytkick Max autosampler (ThermoFisher). Flow cytometry for in vivo experiments Single-cell suspensions of harvested organs were stained with human anti-CD3 BV711 (clone: OKT3), anti-CD4 PerCP-Cy5.5 (clone: OKT4), anti-CD8 PE (clone: HIT8a), anti-CD19 BV786 (clone: HIB19) and anti-CD69 PE-Cy7 (clone: FN50) or mouse anti- CD45 pacific blue (clone: 30-F11) antibodies (Biolegend, USA). Fixable viability dye (eFluor™ 780, eBioscience, USA) was used to exclude dead cells.
	TISSUE PREPARATION - mechanical disintegration - only for tumor tissue: collagense/DNase digestion (37°C, 30min) - pass through a cell strainer to get single cell suspensions - only for brain tissue: brain gradient centrifugation - wash step with PBS/2%FCS - FACS staining (4°C, 30min) - wash step with PBS/2%FCS - FACS analysis in PBS/2% FCS CELL CULTURE EXPERIMENTS - wash T cells with PBS
	- FACS starning (4°C, 30min) - wash with PBS - FACS analysis
Instrument	BD FACS Aria II (BD bioscience, Germany) - for FACS sorting FACS Canto II (BD bioscience, Germany) FACS Fortessa (BD bioscience, Germany) Beckmann Coulter CytoFLEX Thermo Fisher Attune Nxt with autoloader.
Software	Flow cytometric data were analyzed using FlowJo V10.3 to V10.8.1 software, BD FACSDiva (BD bioscience, Germany), and with Matlab 2021b. Data was then processed with Matlab 2021b, Excel 365 or GraphPad Prism V.9.4.0.
Cell population abundance	Post-sort purity was >90%

Binding : FSC/SSC -> Singlet -> Cy5+ signal -> Binding of PTE

Killing: CSFE+ cells -> target cells FSC/SSC -> Dead/alive cells

Activation CSFE- cells -> PBMCs FSC/SCC -> T-cells CD8+ -> CD69+ CD4+ -> CD69+

-> CD8+ -> CD69+

Animal experiments: FSC/SSC -> Singlet -> FVD -> murine CD45 negative -> GFP+ = Nalm-6 GFP+ tumor cells --> Cy5+ : Binding of PTE -> CD19+: Binding of PTE and Blina-BS -> CD3+ -> T cells -> CD4+ -> CD69+

Tick this box to confirm that a figure exemplifying the gating strategy is provided in the Supplementary Information.